



MIC Test Strip Technical Sheet **Mould**

Specimen

Blood, CSF, sterile body fluids and tissues, nasopharynx, respiratory, urinary and gastrointestinal tracts.

| | <i>Aspergillus</i> spp. | <i>Fusarium</i> spp. | <i>Rhizopus</i> spp. | other species |
|----------------------------------|---|--|--------------------------|---|
| Medium | RPMI Agar (Ref. 11509) | | | |
| Inoculum | Suspension of conidia and hyphae (mature growth 5-7 days) in saline with Tween 20 (Ref. 80032) | | | |
| | 0.5 McFarland (Ref. 80400) | 1 McFarland (Ref. 80401) | 1 McFarland (Ref. 80401) | |
| Incubation | 35 ± 2°C / moist ambient air (plates in plastic bag) / 16-72 hours, depending on the genus. | | | |
| | 16-24 hours. | 35 ± 2°C / 24-48 hours followed by room temperature for another 24-48 hours. | 16-24 hours. | Extend incubation time as needed and inspect plates daily for growth and presence of an inhibition ellipse. |
| Interpretation of results | Interpret results after 16-24 hours, confirm at 48 hours if needed. Slow growers may need up to 72 hours. Read the M.I.C. where the inhibition ellipse intersects the M.I.C. scale. For amphotericin B, read complete inhibition of growth and azoles at 80% inhibition. Ignore filaments bending into the ellipse which may be caused by overgrowth in case of prolonged incubation. | | | |

| | Quality Control (MIC µg/mL) 48 hours incubation | | | CLSI INTERPRETATION MIC Criteria (µg/mL) | | | EUCAST INTERPRETATION MIC Criteria (µg/mL) | | Example of ANTIBIOGRAM |
|---|---|----------------------------------|--|---|---|----|---|----------------------|---------------------------|
| | <i>C. parapsilosis</i> ATCC® 22019 | <i>A. flavus</i> ATCC® 204304 | <i>A. fumigatus</i> ATCC® MY A-3626 | S | I | R | S | R | |
| AMB AMPHOTERICIN B <i>A. fumigatus</i> <i>A. niger</i> | 0.25-1 | 0.5-4 | 0.5-4 | | | | ≤1 ≤1 | >2 >2 | ✓ |
| FLU FLUCONAZOLE | 1-8 | - | - | | | | | | |
| FC FLUCYTOSINE | 0.064-0.25 | - | - | | | | | | |
| ITC ITRACONAZOLE <i>A. flavus</i> <i>A. fumigatus</i> <i>A. nidulans</i> <i>A. terreus</i> | 0.064-0.25 | 0.25-0.5 | 0.25-2 | | | | ≤1 ≤1 ≤1 ≤1 | >2 >2 >2 >2 | ✓ |
| KE KETOCONAZOLE | 0.032-0.125 | - | - | | | | | | |
| POS POSACONAZOLE <i>A. fumigatus</i> <i>A. terreus</i> | 0.032-0.25 | 0.06-0.5 | - | | | | ≤0.12 ≤0.12 | >0.25 >0.25 | |
| VO VORICONAZOLE <i>A. fumigatus</i> | 0.016-0.064 | 0.5-4 | 0.25-1 | ≤0.5 | 1 | ≥2 | ≤1 | >2 | |

Susceptible (S), Intermediate (I), Resistant (R)

Note: Antifungal MTS quality control data are not identical to CLSI or EUCAST specifications in all cases. MTS ranges at 48 hours are based on extensive data generated from in house testing.

Disclaimer: The table is intended for general guidance only and may not contain all the necessary information. Also reported interpretive criteria and QC MIC ranges might be out of date. Always current guidelines from CLSI and/or EUCAST should be consulted.

References

- CLSI M61. Performance Standards for Antifungal Susceptibility Testing of Filamentous Fungi. 2nd Edition, 2020.
- EUCAST. Antifungal Agents. Breakpoint tables for interpretation of MICs. Version 9.0, 2018.

CLSI is a trademark belonging to Clinical Laboratory and Standards Institute, Inc.

The ATCC trademark and trade name and any and all ATCC catalog numbers are trademarks of the American Type Culture Collection.

This document has been produced in part under ECDC service contracts and made available at no cost by EUCAST at no cost to the user and can be accessed on the EUCAST website: www.eucast.org. EUCAST recommendations are frequently updated and the latest versions are available at www.eucast.org.

Any other name or trademark is the property of its respective owner.

MTS™ (MIC Test Strip), International Patent

Liofilchem®, the Liofilchem company logo and MTS logo are registered trademarks of LIOFILCHEM s.r.l.



LIOFILCHEM® s.r.l.

Via Scozia, 64026 Roseto degli Abruzzi (TE) Italy
Tel. +39 0858930745 Fax +39 0858930330 www.liofilchem.com

