



# MIC Test Strip Technical Sheet Yeast

## *Candida* spp. and *Cryptococcus neoformans*

### Specimen

Blood, CSF, sterile body fluids and tissues, nasopharynx, urinary, respiratory and gastrointestinal tract.

|                                  | <i>Candida</i> spp.   | <i>Cryptococcus neoformans</i>                                     |
|----------------------------------|---|--|
| <b>Medium</b>                    | RPMI Agar (Ref.11509)   |  |
| <b>Inoculum</b>                  | Suspension in physiological solution  |  |
|                                  | 0.5 McFarland (Ref. 80400)  | 1 McFarland (Ref. 80401)   |
|                                  | Ensure to "double-dip" when inoculating plates i.e. after swabbing the plate the first time, soak the swab and streak the plate again.  |  |
| <b>Incubation</b>                | 35 ± 2°C / moist ambient air (plates in plastic bag) / 24-48 hours.<br>Confirm <i>C. glabrata</i> and <i>C. tropicalis</i> after 48 hours.  | 35 ± 2°C / moist ambient air (plates in plastic bag) / 48-72 hours |
| <b>Interpretation of results</b> | Amphotericin B: interpret at complete inhibition of all growth.<br>Flucytosine: interpret at almost complete (90%) inhibition.<br>Azoles: interpret at the first point of significant inhibition/marked decrease in growth density. Refer to 80% inhibition principle to visually select the M.I.C. endpoint.<br>Echinocandins: interpret trailing endpoints at the first visual point of significant inhibition i.e. 80% inhibition. |  |

|  | Quality Control<br>(MIC µg/mL)<br>48 hours incubation |                                |                                   | CLSI INTERPRETATION<br>MIC Criteria (µg/mL)  |                                     |                                     | EUCAST INTERPRETATION<br>MIC Criteria (µg/mL) |  | Example of<br>ANTIBIOGRAM<br><br>140 mm petri dish |
|--|---|--------------------------------|-----------------------------------|--|-------------------------------------|-------------------------------------|---|--|--|
|  | <i>C. parapsilosis</i><br>ATCC® 22019                 | <i>C. krusei</i><br>ATCC® 6258 | <i>C. albicans</i><br>ATCC® 90028 | S  | I                                   | R                                   | S   | R                                      |  |
| <b>AMB</b> AMPHOTERICIN B<br><i>C. albicans</i><br><i>C. glabrata</i><br><i>C. tropicalis</i><br><i>C. krusei</i><br><i>C. parapsilosis</i>                            | 0.25-1  | 0.5-2                          | 0.125-0.5                         |  |                                     |                                     | ≤1<br>≤1<br>≤1<br>≤1<br>≤1                    | >1<br>>1<br>>1<br>>1<br>>1             | ✓  |
| <b>AND</b> ANIDULAFUNGIN<br><i>C. albicans</i><br><i>C. glabrata</i><br><i>C. tropicalis</i><br><i>C. krusei</i><br><i>C. parapsilosis</i><br><i>C. guilliermondii</i> | 0.5-4   | 0.016-0.125                    | 0.002-0.008                       | ≤0.25<br>≤0.12<br>≤0.25<br>≤0.25<br>≤2<br>≤2 | 0.5<br>0.25<br>0.5<br>0.5<br>4<br>4 | ≥1<br>≥0.5<br>≥1<br>≥1<br>≥8<br>≥8  | ≤0.03<br>≤0.06<br>≤0.06<br>≤0.06<br>≤0.002    | >0.03<br>>0.06<br>>0.06<br>>0.06<br>>4 |  |
| <b>CAS</b> CASPOFUNGIN<br><i>C. albicans</i><br><i>C. glabrata</i><br><i>C. tropicalis</i><br><i>C. krusei</i><br><i>C. parapsilosis</i><br><i>C. guilliermondii</i>   | 0.25-2  | 0.25-1                         | 0.064-0.25                        | ≤0.25<br>≤0.12<br>≤0.25<br>≤0.25<br>≤2<br>≤2 | 0.5<br>0.25<br>0.5<br>0.5<br>4<br>4 | ≥1<br>≥0.5<br>≥1<br>≥1<br>≥8<br>≥8  |   |  |  |
| <b>FLU</b> FLUCONAZOLE<br><i>C. albicans</i><br><i>C. glabrata</i><br><i>C. parapsilosis</i><br><i>C. tropicalis</i>   | 1-8   | 128-≥256                       | 0.125-0.5                         | ≤2<br>-<br>≤2<br>≤2                          | 4<br>32<br>4<br>4                   | ≥8<br>≥64<br>≥8<br>≥8               | ≤2<br>≤0.002<br>≤2<br>≤2                      | >4<br>>32<br>>4<br>>4                  | ✓  |
| <b>FC</b> FLUCYTOSINE<br><i>C. albicans</i><br><i>C. dubliniensis</i><br><i>C. parapsilosis</i><br><i>C. tropicalis</i>  | 0.064-0.25  | ≥32                            | 0.5-2                             |  |                                     |                                     |   |  | ✓  |
| <b>ITC</b> ITRACONAZOLE<br><i>C. albicans</i><br><i>C. dubliniensis</i><br><i>C. parapsilosis</i><br><i>C. tropicalis</i>  | 0.064-0.25  | 0.25-1                         | 0.064-0.25                        |  |                                     |                                     | ≤0.06<br>≤0.06<br>≤0.12<br>≤0.12              | >0.06<br>>0.06<br>>0.12<br>>0.12       | ✓  |
| <b>KE</b> KETOCONAZOLE<br><i>C. albicans</i><br><i>C. dubliniensis</i><br><i>C. parapsilosis</i><br><i>C. tropicalis</i>   | 0.032-0.125   | 0.25-1                         | 0.008-0.032                       |  |                                     |                                     |   |  |  |
| <b>MYC</b> MICAFAUNGIN<br><i>C. albicans</i><br><i>C. glabrata</i><br><i>C. tropicalis</i><br><i>C. krusei</i><br><i>C. parapsilosis</i><br><i>C. guilliermondii</i>   | 0.25-2  | 0.032-0.25                     | 0.004-0.032                       | ≤0.25<br>≤0.06<br>≤0.25<br>≤0.25<br>≤2<br>≤2 | 0.5<br>0.12<br>0.5<br>0.5<br>4<br>4 | ≥1<br>≥0.25<br>≥1<br>≥1<br>≥8<br>≥8 | ≤0.016<br>≤0.03<br>≤0.002                     | >0.016<br>>0.03<br>>2                  |  |
| <b>POS</b> POSACONAZOLE<br><i>C. albicans</i><br><i>C. dubliniensis</i><br><i>C. parapsilosis</i><br><i>C. tropicalis</i>  | 0.032-0.25  | 0.125-0.5                      | 0.032-0.125                       |  |                                     |                                     | ≤0.06<br>≤0.06<br>≤0.06<br>≤0.06              | >0.06<br>>0.06<br>>0.06<br>>0.06       |  |

|                        | Quality Control<br>(MIC µg/mL)<br>48 hours incubation | CLSI INTERPRETATION<br>MIC Criteria (µg/mL) |                                   |       | EUCAST INTERPRETATION<br>MIC Criteria (µg/mL) |    | Example of<br>ANTIBIOGRAM |       |
|------------------------|---|---|-----------------------------------|-------|---|----|---------------------------|-------|
|                        |   | S   | I                                 | R     | S   | R  |                           |       |
| <b>VO</b> VORICONAZOLE | <i>C. parapsilosis</i><br>ATCC® 22019                 | <i>C. krusei</i><br>ATCC® 6258              | <i>C. albicans</i><br>ATCC® 90028 |       |   |    | 140 mm petri dish         |       |
|                        | 0.016-0.064   | 0.25-1                                      | 0.004-0.016                       |       |   |    | ✓                         |       |
| <i>C. albicans</i>     |   |   |                                   | ≤0.12 | 0.25-0.5                                      | ≥1 | ≤0.06                     | >0.25 |
| <i>C. krusei</i>       |   |   |                                   | 0.5   | 1   | ≥2 | ≤0.06                     | >0.25 |
| <i>C. parapsilosis</i> |   |   |                                   | ≤0.12 | 0.25-0.5                                      | ≥1 | ≤0.12                     | >0.25 |
| <i>C. tropicalis</i>   |   |   |                                   | ≤0.12 | 0.25-0.5                                      | ≥1 | ≤0.12                     | >0.25 |

Susceptible (S), Intermediate (I), Resistant (R)

**Note:** Antifungal MTS quality control data are not identical to CLSI or EUCAST specifications in all cases. MTS ranges at 48 hours are based on extensive data generated from in house testing.

**Disclaimer:** The table is intended for general guidance only and may not contain all the necessary information. Also reported interpretive criteria and QC MIC ranges might be out of date. Always current guidelines from CLSI and/or EUCAST should be consulted.

#### References

1. CLSI M60. Performance Standards for Antifungal Susceptibility Testing of Yeasts. 2nd Edition, 2020.
2. EUCAST. Antifungal Agents. Breakpoint tables for interpretation of MICs. Version 9.0, 2018.

CLSI is a trademark belonging to Clinical Laboratory and Standards Institute, Inc.

The ATCC trademark and trade name and any and all ATCC catalog numbers are trademarks of the American Type Culture Collection.

This document has been produced in part under ECDC service contracts and made available at no cost by EUCAST at no cost to the user and can be accessed on the EUCAST website: [www.eucast.org](http://www.eucast.org). EUCAST recommendations are frequently updated and the latest versions are available at [www.eucast.org](http://www.eucast.org).

Any other name or trademark is the property of its respective owner.

#### MTS™ (MIC Test Strip), International Patent

Liofilchem®, the Liofilchem company logo and MTS logo are registered trademarks of LIOFILCHEM s.r.l.



**LIOFILCHEM® s.r.l.**

Via Scozia, 64026 Roseto degli Abruzzi (TE) Italy  
Tel. +39 0858930745 Fax +39 0858930330 [www.liofilchem.com](http://www.liofilchem.com)

