MTS™

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Quantitative assay for determining the Minimum Inhibitory Concentration (MIC).

DESCRIPTION

MTS™ is a quantitative assay for determining the Minimum Inhibitory Concentration (MIC) of antimicrobial agents against microorganisms and for detecting the resistance mechanisms. **MTS™** are paper strips with special features* that are impregnated with a predefined concentration gradient of antibiotic, across 15 two-fold dilutions of a conventional MIC method.

On one side of the strip is indicated a MIC scale in µg/mL and a code that identify the antimicrobial agent.

For ESBL (Extended Spectrum Beta Lactamase) and MBL (Metallo Beta Lactamase) detection, the double-sided gradient carries the appropriate diagnostic reagents.

MTS™ are available in a large variety of configurations. Each configuration is available in packages of 10, 30 and 100 tests.

CONTENTS OF THE PACKAGES

The 10-test box contains 10 strips individually packed in desiccant envelopes and an instruction sheet.

The 30-test box contains 30 strips individually packed in desiccant envelopes and an instruction sheet.

The 100-test box contains 10 desiccant envelopes, each containing 10 strips, and an instruction sheet. The 100-test pack also contains a storage tube.

METHOD PRINCIPLE

When the **MTS™** is applied onto an inoculated agar surface, the preformed exponential gradient of antimicrobial agent is transferred to the agar matrix.

After 18 hours incubation or longer, a symmetrical inhibition ellipse centered along the strip is formed. The MIC is read directly from the scale in terms of µg/mL at the point where the edge of the inhibition ellipse intersects the strip **MTS™**.

Other growth/inhibition patterns may also be seen for resistance detection methods.

COMPOSITION

The strips are made of high-quality paper and each strip is impregnated with a predefined concentration gradient across 15 two-fold dilutions of antibiotic agent.

GATHERING AND KEEPING SAMPLES

The colonies that are to be subjected to the evaluation of Minimum Inhibition Concentration (MIC) are taken up by culture media that have been previously swabbed with the sample under examination. In the case of mixed colonies the bacterial strains must be purified before inoculation.

TEST PROCEDURE

1. Allow unopened envelop to come to room temperature before opening it, for minimising condensation on the strip.
2. Swab 4 to 5 well isolated and morphologically similar colonies with a culture medium and suspend them in 5 mL of a suitable suspension medium. Fastidious microorganisms should be suspended in broth and used within 15 minutes.
3. Compare the turbidity to the appropriate McFarland standard.
4. Dip a sterile swab in the broth culture or in a diluted form thereof and squeeze it on the wall of the test tube to eliminate excess liquid.
5. Drag it along the surface of the medium contained on the plate so as to produce even growth; allow excess moisture to be absorbed and ensure that the surface is completely dry before applying strips.
6. Apply the strip to the agar surface with the MIC scale facing upwards and code of the strip to the outside of the plate, pressing it with a sterile forceps on the surface of the agar and ensure that whole length of the antibiotic gradient is in complete contact with the agar surface. Once applied, do not move the strip.
7. Incubate plates in an inverted position under conditions appropriate for the microorganism.
8. Put the not used strips onto the tube contained in the package.

EVALUATING THE RESULTS

At the end of incubation read the MIC value where the edge of the inhibition ellipse intersects the strip (intersection between two scale segments should be round up to the higher value).

MIC break points for defining susceptibility categories as provided by CLSI or EUCAST could be used for interpreting MIC values.

Always round up **MTS™** half dilution values to the next upper two-fold value before categorisation. An overview of CLSI and EUCAST interpretative criteria is provided in Table no.1.

CLINICAL INTERPRETATION

The test **MTS™** carried out *in vitro* cannot exactly reproduce *in vivo* conditions. Nevertheless, it shows the effect of the concentration of the antibiotic, which varies in the culture medium in relation to the growth of the microbial population. The final choice of antibiotic to administer to the patient is the responsibility of the clinician who possesses all the information on the patient.

QUALITY CONTROL

Each batch of **MTS™** is subjected to precise and thorough checks in compliance with CLSI standards using the bacterial strains indicated in the table no.1.

PRECAUTIONS

The **MTS™** cannot be classified as being hazardous according to current legislation. **MTS™** are disposable products. **MTS™** are only for diagnostic *in vitro* use and are intended for professional use. They must be used in the laboratory by properly trained operators using approved aseptic and safety methods for pathogenic agents.

STORAGE

The unopened package of **MTS™** should be stored at -20°C until the given expiry date.

Leftover **MTS™** from an opened package must be stored at 2-8°C in the airtight tube, containing desiccant, provided in the pack for no more than 7 days.

Do not store near sources of heat and do not expose to excessive temperature variations.

ELIMINATING USED MATERIAL

After use, **MTS™** and the material that comes into contact with the sample must be decontaminated and disposed of in accordance with current laboratory techniques for the decontamination and disposal of potentially infected material.

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- CLSI. Performance Standards for Antimicrobial Susceptibility Testing; 29th ed. CLSI Supplement M100S. Wayne, PA: Clinical and Laboratory Standards Institute; 2019.
 - CLSI. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically; 11th ed. CLSI standard M07. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
 - CLSI. Methods for Dilution Antimicrobial Susceptibility Testing of Anaerobic Bacteria; 9th ed. CLSI standard M11. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
 - CLSI. Performance Standards for Antifungal Susceptibility Testing of Yeasts; 1st ed. CLSI Supplement M60. Wayne, PA: Clinical and Laboratory Standards Institute; 2017.
 - CLSI. Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts; 4th ed. CLSI standard M27. Wayne, PA: Clinical and Laboratory Standards Institute; 2017.
 - CLSI. Performance Standards for Antifungal Susceptibility Testing of Filamentous Fungi; 1st ed. CLSI Supplement M61. Wayne, PA: Clinical and Laboratory Standards Institute; 2017.
 - CLSI. Reference Method for Broth Dilution Antifungal Susceptibility Testing of Filamentous Fungi; 3rd ed. CLSI standard M38. Wayne, PA: Clinical and Laboratory Standards Institute; 2017.
 - CLSI. Performance Standards for Antimicrobial Disk and Dilution Susceptibility Tests for Bacteria Isolated From Animals; 5th ed. CLSI standard VET01. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
 - CLSI. Performance Standards for Antimicrobial Disk and Dilution Susceptibility Tests for Bacteria Isolated From Animals; Approved Standard - Fourth Edition. CLSI document VET01-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
 - CLSI. Performance Standards for Antimicrobial Susceptibility Testing of Bacteria Isolated From Aquatic Animals; Second Informational. CLSI document VET03/VET04-S2. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
 - CLSI. Methods for Broth Dilution Susceptibility Testing of Bacteria Isolated From Aquatic Animals; Approved Guideline - Second Edition. CLSI document VET04-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
 - CLSI. Methods for Antimicrobial Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria Isolated From Animals; 1st ed. CLSI supplement VET06. Wayne, PA: Clinical and Laboratory Standards Institute; 2017.
 - The European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MICs and zone diameters. Version 9.0, 2019. <http://www.eucast.org>
 - The European Committee on Antimicrobial Susceptibility Testing. Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST. Version 9.0, 2019. <http://www.eucast.org>
 - The European Committee on Antimicrobial Susceptibility Testing. Antifungal Agents. Breakpoint tables for interpretation of MICs. Version 9.0, 2018. <http://www.eucast.org>
-
- Widlake, E. et al. (2019). Evaluation of three methodologies for in vitro susceptibility testing of Ceftolozane-Tazobacatam (C/T). FEMS Microbe poster 183.
 - Olesky M. et al. (2019). Multi-Site Evaluation of Eravacycline MIC Test Strip (MTS™) Compared To Broth Microdilution MICs. MAD-ID poster.
 - Koeth, L. M. et al. (2018). Comparison of Plazomicin MIC Test Strip and Broth Microdilution MIC Results for 125 Enterobacteriaceae. ID Week poster 2060.
 - Koeth, L. M. et al. (2018). Multi-Site Evaluation of Meropenem/Vaborbactam MIC Test Strip (MTS) Compared To Broth Microdilution MICs. ASM Microbe poster 240.
 - Savini, V. et al. (2016). Daptomycin-resistant *Staphylococcus pettenkoferi* of human origin. Acta Biochim Pol; 63(2): 297-301.
 - Koeth, L. M. et al. (2016). Evaluation of Ceftolozane-Tazobactam MTS™ Compared to Broth Microdilution MIC for Enterobacteriaceae and *Pseudomonas aeruginosa*. ASM Microbe poster.
 - Koeth, L. M. et al. (2016). Multi-Site Evaluation of Dalbavancin and Vancomycin MTS™ Compared To Broth Microdilution MICs. ASM Microbe poster.
 - Carreto, E. et al. (2016). Multicentric evaluation of the reliability and the reproducibility of synergy testing using the MTS™ - synergy application system (MTS-SAS™). ECCMID poster 1009.
 - Koeth, L. M. et al. (2015) Evaluation of Dalbavancin MTS™ (MTS) Compared to Broth Microdilution MIC for Relevant Gram Positive Isolates. ICAAC poster D-1138.
 - Savini, V. et al. (2014). Liofilchem® Chromatic VRE and vancomycin MTS™ detected glycopeptide resistance in a vanB neonatal *Enterococcus faecium* isolate showing alternate vancomycin susceptibility and resistance with bioMérieux Vitek2. Int J Clin Exp Pathol; 7(9):6274-6277.
 - Canton, R. et al. (2014). Determination of the sensitivity and specificity of the MTS™s (MTS) mechanisms of detection against in-house methods for 644 multi-drug resistant (MDR) strains: a European multi-centre study. ECCMID poster 1009.

- D'Humières, C. et al (2014). Direct determination of antimicrobial susceptibility of Gram Negative Bacilli from respiratory samples on chromogenic agar plates with gradient antibiotic strips. ECCMID poster 348.
- Tyrrell, J. and Walsh, T.R. (2014). Efficacy of Fosmidomycin alone and in combination with colistin, tigecycline and rifampicin against multi-drug resistant and extensively drug resistant Enterobacteriaceae. ECCMID poster 235.
- Walsh, T. et al. (2013). Synergistic effects of rifampicin, nitrofurantoin, fosfomycin and colistin against NDM-1 positive *Klebsiella pneumoniae* and *Escherichia coli* using the new MTS synergy method. ECCMID poster 1553.
- Walsh, T.R. et al. (2013). Evaluation of MTSTM's (MTS) gradient-diffusion system for susceptibility testing of NDM-1 positive *Klebsiella pneumoniae* and *Escherichia coli*. ECCMID poster 1554.

- Savini, V. et al. (2013). Methicillin-Resistant *Staphylococcus pseudintermedius* Infection in a Bone Marrow Transplant Recipient. *J Clin Microbiol*; 51(5):1636.
- Meletiadis, J. et al. (2013). Evaluation of gradient concentration strips for in vitro combination testing of antifungal combinations against *Candida* spp. ECCMID poster 976.
- Meletiadis, J. et al. (2013). Evaluation of strips with concentration gradient of antifungal agents for the susceptibility determination of *Candida* strains ECCMID poster 1583.
- Brocco, F. et al. (2013). Evaluation of a new method for detecting KPC-producing bacteria based on the combination of a screening medium and a gradient diffusion system. ECCMID poster eP690.
- M'Zali, F.H. et al. (2013). A novel, direct susceptibility testing method of sonicated vascular prosthetic graft samples by combination of Liofilchem® MTSTM and Chromatic™ agar plates. ECCMID poster 1562.
- Haldorsen, B. C. et al. (2012). Agreement of the MTSTM vs. Etest in the MIC determination of *Streptococcus pneumoniae*. ECCMID poster 1658.
- Rossolini, G.M. et al. (2011). Evaluation of a new gradient-diffusion system for MIC determination with Gram-negative pathogens. ECCMID, poster 572.
- Stefani, S. et al (2011). A new reliable screening method for the evaluation of VISA and hVISA strains by "Vancomycin-Teicoplanin MTS™" (VTMTS). ECCMID poster 776.

TABELLA DEI SIMBOLI / TABLE OF SYMBOLS / TABLEAU DES SYMBOLES / SYMBOLE / TABLA DE SIMBOLOS / TABELA DE SÍMBOLOS / SYMBOLFORKLARINGER / TABUĽKA SYMBOLOV / TABULKA SYMBOLŮ / SYMBOLTAVLE / ΠΙΝΑΚΑΣ ΣΥΜΒΟΛΩΝ / SYMBOLER

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Liofilchem.com/MTS

MTS™ (MIC Test Strip), International Patent

Liofilchem®, the Liofilchem company logo and the MTS logo are registered trademarks of LIOFILCHEM s.r.l.

© Copyright LIOFILCHEM 2020

LIOFILCHEM® s.r.l.

Via Scozia zona ind.le, 64026 Roseto degli Abruzzi (Te) Italy
Tel. +39 0858930745 Fax +39 0858930330 www.liofilchem.com

In USA, available for products noted as "FDA Cleared" in the MTS™ Catalog.

F00023
Rev.35 / 03.01.2020