# EUCAST ESBL disc kit

ENGLISH

Disc tests for confirmation of ESBL in Enterobacterales.

## DESCRIPTION

Extended-spectrum  $\beta$ -lactamases (ESBLs) are enzymes hydrolyzing most penicillins and cephalosporins, including oxyimino- $\beta$ -lactam compounds but not cephamycins and carbapenems. Most ESBLs belong to the Ambler class A of  $\beta$ -lactamases and are inhibited by  $\beta$ -lactamases inhibitors: clavulanic acid, sulbactam and tazobactam. ESBL production has been observed mostly in Enterobacterales, particularly *Escherichia coli* and *Klebsiella pneumoniae*, but all other clinically-relevant Enterobacterales species are also common ESBL-producers. In many areas, ESBL detection and characterization is recommended or mandatory for infection control purpose.

ESBL detection involves two important steps. The first is a screening test with an indicator cephalosporin which looks for resistance or diminished susceptibility, thus identifying isolates likely to be harboring ESBLs. The second step is a confirmation test which evaluates the synergy between an oxyimino cephalosporin and clavulanic acid, distinguishing isolates with ESBLs from those that are resistant for other reasons.

## CONTENTS OF THE PACKAGES

6 x 50 discs cartridges, each packaged in a "blister" with a dryer.

## METHOD PRINCIPLE

Enterobacterales suspected to be producers of ESBL enzymes may be confirmed by evaluating the inhibition of ESBL activity by **Clavulanic acid**.

### Combination Disc Test (CDT)

For each test, discs containing cephalosporin alone (cefotaxime AND ceftazidime, OR cefepime) and in combination with clavulanic acid are applied. The inhibition zone around the cephalosporin disc combined with clavulanic acid is compared with the zone around the disc with the cephalosporin alone. The test is positive if the inhibition zone diameter is  $\geq 5$  mm larger with clavulanic acid than without.

When high level expression of AmpC  $\beta$ -lactamases is suspected, use cefepime as the indicator cephalosporin, as cefepime is not hydrolyzed by AmpC  $\beta$ -lactamases. Resistance to cephamycins, e.g. a ceftaxime MIC  $> 8$   $\mu\text{g/ml}$ , may be indicative of high level expression of AmpC  $\beta$ -lactamases.

## GATHERING AND KEEPING SAMPLES

The colonies that are to be subjected to the susceptibility test are taken up by culture media that have been previously swabbed with the sample under examination.

## TEST PROCEDURE

- Using a fresh, pure culture prepare a suspension of the test organism equal to 0.5 McFarland Standard.
- Using a sterile cotton swab, spread the adjusted suspension over the entire area of a Mueller Hinton agar plate.
- Apply the discs onto the inoculated plate, ensuring sufficient space between individual discs to allow for proper measurement of inhibition zones.
- Incubate at  $35 \pm 1^\circ\text{C}$  for 16-20 hours.

## EVALUATING THE RESULTS

At the end of the incubation period, measure the inhibition halos and interpret as indicated in the following table.

### Interpretative Table. CDT method for confirmation of ESBL.

Antibiotic Disc	Confirmation is positive if
Cefotaxime alone and with Clavulanic acid *    CTX 5 $\mu\text{g}$ and CTL 5+10 $\mu\text{g}$ AND/OR	$\geq 5$ mm increase in inhibition zone of cephalosporin with Clavulanic acid
Ceftazidime alone and with Clavulanic acid    CAZ 10 $\mu\text{g}$ and CAL 10+10 $\mu\text{g}$	
Cefepime alone and with Clavulanic acid **    FEP 30 $\mu\text{g}$ and FEL 30+10 $\mu\text{g}$	
* <i>Escherichia coli</i> , <i>Klebsiella</i> spp., <i>Proteus mirabilis</i> , <i>Salmonella</i> spp., <i>Shigella</i> spp.	
** <i>Enterobacter</i> spp., <i>Citrobacter freundii</i> , <i>Morganella morganii</i> , <i>Providencia</i> spp., <i>Hafnia alvei</i> (all producing chromosomal AmpC)	

## QUALITY CONTROL

To check the performance of discs, MHA plates, inoculum and procedure used, test *Escherichia coli* ATCC® 25922 as a negative control for ESBL. As a positive control either *Klebsiella pneumoniae* ATCC® 700603 (SHV-18 ESBL) or *Escherichia coli* NCTC 13353 (CTX-M-15 ESBL) can be used, depending on the combination disc. See tables on the following page for specific criteria. Target values and ranges for Cefotaxime 5  $\mu\text{g}$ , Ceftazidime 10  $\mu\text{g}$  and Cefepime 30  $\mu\text{g}$  can be found in the EUCAST QC Tables available at [www.eucast.org](http://www.eucast.org).

Control strain	Cefotaxime (5 µg)	Ceftazidime (10 µg)	Result
	Increase in zone diameter with Clavulanic acid (10 µg) compared with the cephalosporin alone		
<i>Klebsiella pneumoniae</i> ATCC® 700603 <sup>a</sup>	≥ 3 mm	≥ 5 mm	Positive
<i>Escherichia coli</i> NCTC 13353 <sup>a,b</sup>	≥ 3 mm	≥ 5 mm	Positive
<i>Escherichia coli</i> ATCC® 25922	± 1 mm	± 1 mm	Negative

Control strain	Increase in zone diameter with Cefepime (30 µg) + Clavulanic acid (10 µg) compared with Cefepime alone	Result
<i>Escherichia coli</i> NCTC 13353 <sup>b,c</sup>	≥ 5 mm	Positive
<i>Escherichia coli</i> ATCC® 25922	± 1 mm	Negative

<sup>a</sup> Strain free of choice; one of the strains has to be used as a minimum.

<sup>b</sup> CTX-M-15-producing strain *E. coli* NCTC 13353, which is recommended by CLSI for quality control of various cephalosporins and beta-lactam combination agents, is not included in the EUCAST guidelines.

<sup>c</sup> Zone diameter QC range for Cefepime 30 µg is 6-15 mm; If discrete colonies or a haze of growth are present inside the zone of inhibition, measure the colony-free inner zone (CLSI M100 ED33:2023).

## LIMITS

Diffusion susceptibility tests use an *in vitro* technique and cannot therefore reproduce the extremely complex *in vivo* conditions. Nevertheless, it is a useful and important tool that helps the clinician choose the correct therapy. Many variable factors influence the final result of the diffusion susceptibility test. The main ones are: the culture medium used, impregnation of the discs, inoculation of the medium, temperature, time and incubation atmosphere of the plates, pre-incubation and pre-diffusion conditions, depth of the medium, etc.

## PRECAUTIONS

The disc cannot be classified as being hazardous according to current legislation but fall within the specific field of application where a safety data sheet must be supplied because they can cause phenomena of sensitization in sensitive subjects if they come into contact with the skin. The discs are disposable products for *in-vitro* diagnostic use (\*). They must be used in the laboratory by properly trained operators using approved aseptic and safety methods for pathogenic agents.

\* In USA, these products are available for Research Use Only (RUO) and should not be used for diagnostic purposes.

## STORAGE

Store the unopened blister at -20°C to +8°C till the expiry date. Allow unopened cartridge to come to room temperature before removing it from the blister for minimizing condensation on the discs. Leftover discs from an opened cartridge should be stored at 2-8°C for no more than 7 days. Return unused discs to the refrigerator as soon as the application of the discs has been completed. Dispose of expire discs.

## ELIMINATING USED MATERIAL


After use, the discs and the material that comes into contact with the sample must be decontaminated and disposed of in accordance with current laboratory techniques for the decontamination and disposal of potentially infected material.

## REFERENCES

- EUCAST technical guidance on the use of the combination disk test (CDT) for confirmation of ESBL in Enterobacterales - New disk potencies for combination disks containing cefotaxime and ceftazidime without and with clavulanic acid (12 February 2019).
- EUCAST guidelines for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance. Version 2.0, 2017.

Disc	µg	Disc	µg	Σ	REF
Cefotaxime	CTX 5	Cefotaxime + Clavulanic acid	CTL 15 (5+10)	50 Test	99002
Ceftazidime	CAZ 10	Ceftazidime + Clavulanic acid	CAL 20 (10+10)		
Cefepime	FEP 30	Cefepime + Clavulanic acid	FEL 40 (30+10)		

## TABLE OF SYMBOLS

<b>LOT</b> Batch code	<b>IVD</b> In Vitro Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
<b>REF</b> Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Consult instructions for use	 Do not reuse

## Revision History

Revision	Release Date	Change Summary
1	2024-02-06	Revised QUALITY CONTROL

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**Liofilchem, Inc.** In US available as RUO device