ESBL disc kit (acc. to CLSI)

Disc tests for confirmation of ESBL-producing Enterobacteriaceae.

DESCRIPTION

Extended-spectrum β -lactamases (ESBLs) are enzymes hydrolyzing most penicillins and cephalosporins, including oxyimino- β -lactam compounds but not cephamycins and carbapenems. Most ESBLs belong to the Ambler class A of β -lactamases and are inhibited by β -lactamases inhibitors: clavulanic acid, sulbatam and tazobactam. ESBL production has been observed mostly in Enterobacteriaceae, particularly *Escherichia coli* and *Klebsiella pneumoniae*, but all other clinically-relevant Enterobacteriaceae species are also common ESBL-producers. In many areas, ESBL detection and characterization is recommended or mandatory for infection control purpose.

ESBL detection involves two important steps. The first is a screening test with an indicator cephalosporin which looks for resistance or diminished susceptibility, thus identifying isolates likely to be harboring ESBLs. The second step is a confirmation test which evaluates the synergy between an oxyimino cephalosporin and clavulanic acid, distinguishing isolates with ESBLs from those that are resistant for other reasons.

CONTENTS OF THE PACKAGES

4 x 50 discs cartridges, each packaged in a "blister" with a dryer.

METHOD PRINCIPLE

Enterobacteriaceae suspected to be producers of ESBL enzymes may be confirmed by evaluating the inhibition of ESBL activity by **Clavulanic acid**.

Combination Disc Test (CDT)

For each test, discs containing both cefotaxime and ceftazidime, alone and in combination with clavulanic acid are applied. The inhibition zone around the cephalosporin disc combined with clavulanic acid is compared with the zone around the disc with the cephalosporin alone. The test is positive if the inhibition zone diameter of either cephalosporin is ≥ 5 mm larger with clavulanic acid than without.

GATHERING AND KEEPING SAMPLES

The colonies that are to be subjected to the susceptibility test are taken up by culture media that have been previously swabbed with the sample under examination.

TEST PROCEDURE

1. Using a fresh, pure culture prepare a suspension of the test organism equal to 0.5 McFarland Standard.
2. Using a sterile cotton swab, spread the adjusted suspension over the entire area of a Mueller Hinton agar plate.
3. Apply the discs onto the inoculated plate, ensuring sufficient space between individual discs to allow for proper measurement of inhibition zones.
4. Incubate at $35 \pm 2^\circ\text{C}$ for 18-24 hours.

EVALUATING THE RESULTS

At the end of the incubation period, measure the inhibition halos and interpret as indicated in the following table.

Interpretative Table. CDT method for confirmation of ESBL in *K. pneumoniae*, *K. oxytoca*, *E. coli* and *P. mirabilis*.

Antibiotic Disc		Confirmation is positive if
Cefotaxime alone and with Clavulanic acid	CTX 30 μg and CTL 30+10 μg	≥ 5 mm increase in inhibition zone of either cephalosporin with Clavulanic acid
Ceftazidime alone and with Clavulanic acid	CAZ 30 μg and CAL 30+10 μg	

QUALITY CONTROL

Appropriate strains for quality control of ESBL detection tests:

Microorganism		ESBL phenotype	AmpC phenotype
<i>Klebsiella pneumoniae</i>	ATCC® 700603	+	-
<i>Enterobacater cloacae</i>	ATCC® BAA-1143	-	+
<i>Escherichia coli</i>	ATCC® 25922	-	-

LIMITS

Diffusion susceptibility tests use an *in vitro* technique and cannot therefore reproduce the extremely complex *in vivo* conditions. Nevertheless, it is a useful and important tool that helps the clinician choose the correct therapy. Many variable factors influence the final result of the diffusion susceptibility test. The main ones are: the culture medium used, impregnation of the discs, inoculation of the medium, temperature, time and incubation atmosphere of the plates, pre-incubation and pre-diffusion conditions, depth of the medium, etc.

PRECAUTIONS

The disc cannot be classified as being hazardous according to current legislation but fall within the specific field of application where a safety data sheet must be supplied because they can cause phenomena of sensitization in sensitive subjects if they come into contact with the skin.

The discs are disposable products. They are only for diagnostic *in vitro* use and are intended for professional use. They must be used in the laboratory by properly trained operators using approved aseptic and safety methods for pathogenic agents.

STORAGE

Store the unopened blister at -20°C to +8°C till the expiry date. Allow unopened cartridge to come to room temperature before removing it from the blister for minimising condensation on the discs. Leftover discs from an opened cartridge should be stored at 2-8°C for no more than 7 days. Return unused discs to the refrigerator as soon as the application of the discs has been completed. Dispose of expire discs.

ELIMINATING USED MATERIAL

After use, the discs and the material that comes into contact with the sample must be decontaminated and disposed of in accordance with current laboratory techniques for the decontamination and disposal of potentially infected material.

REFERENCES

- CLSI M100-S24 - Performance Standards for Antimicrobial Susceptibility Testing; 2014.

PRESENTATION

DESCRIPTION		µg	PACKAGING	REF
Cefotaxime	CTX	30	50 Discs	99005
Cefotaxime + Clavulanic acid	CTL	30+10	50 Discs	
Ceftazidime	CAZ	30	50 Discs	
Ceftazidime + Clavulanic acid	CAL	30+10	50 Discs	

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In Vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult accompanying documents	 Do not reuse

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