ESBL (Chromos. Induc. AmpC) disc kit (acc. to EUCAST)

Disc tests for confirmation of ESBLs in organisms with chromosomally encoded inducible AmpC.

DESCRIPTION

Extended-spectrum β -lactamases (ESBLs) are enzymes hydrolyzing most penicillins and cephalosporins, including oxyimino- β -lactam compounds but not cephamycins and carbapenems. Most ESBLs belong to the Ambler class A of β -lactamases and are inhibited by β -lactamase inhibitors: clavulanic acid, sulbactam and tazobactam. ESBL production has been observed mostly in Enterobacteriaceae, particularly *Escherichia coli* and *Klebsiella pneumoniae*, but all other clinically-relevant Enterobacteriaceae species are also common ESBL-producers. In many areas, ESBL detection and characterization is recommended or mandatory for infection control purpose.

AmpC β -lactamases are enzymes encoded on the chromosome of many Enterobacteriaceae and a few other organisms. AmpC enzymes belong to class C, are active on penicillins but even more active on cephalosporins and can significantly hydrolyze cephamycins such as cefoxitin and cefotetan and oxyiminocephalosporins such as ceftazidime, cefotaxime, and ceftriaxone. Inhibitors of class A enzymes such as clavulanic acid, sulbactam, and tazobactam have much less effect on AmpC, although some are inhibited by tazobactam or sulbactam. In many Enterobacteriaceae, AmpC expression is low but inducible in response to β -lactam exposure. β -lactamase inhibitors are also inducers, especially clavulanate. Overexpression confers resistance to broad-spectrum cephalosporins and is a problem especially in infections due to *E. aerogenes* and *E. cloacae*, where an isolate initially susceptible to these agents may become resistant upon therapy. Cefepime is a poor inducer of AmpC β -lactamase, rapidly penetrates through the outer cell membrane, and is little hydrolyzed by the enzyme, so many AmpC-producing organisms test cefepime susceptible.

ESBL detection involves two important steps. The first is a screening test with an indicator cephalosporin which looks for resistance or diminished susceptibility, thus identifying isolates likely to be harboring ESBLs. The second step is a confirmation test which evaluates the synergy between an oxyimino cephalosporin and clavulanic acid, distinguishing isolates with ESBLs from those that are resistant for other reasons. On the contrary, organisms producing enough AmpC β -lactamase will typically give a positive ESBL screening test but fail the confirmatory test involving increased sensitivity with clavulanic acid. To recognize organisms which simultaneously produce both ESBL and AmpC enzymes cefepime is the antibiotic of choice.

CONTENTS OF THE PACKAGES

4 x 50 discs cartridges, each packaged in a "blister" with a dryer.

METHOD PRINCIPLE

Enterobacteriaceae suspected to be producers of both ESBL and AmpC enzymes may be confirmed by using cefepime to evaluate the inhibition of ESBL activity by **Clavulanic acid**.

Combination Disc Test (CDT)

For each test, discs containing cefepime alone and in combination with clavulanic acid are applied. The inhibition zone around the cefepime disc combined with clavulanic acid is compared with the zone around the disc with the cefepime alone. The test is positive if the inhibition zone diameter is ≥ 5 mm larger with clavulanic acid than without.

GATHERING AND KEEPING SAMPLES

The colonies that are to be subjected to the susceptibility test are taken up by culture media that have been previously swabbed with the sample under examination.

TEST PROCEDURE

1. Using a fresh, pure culture prepare a suspension of the test organism equivalent to 0.5 McFarland.
2. Using a sterile cotton swab, spread the adjusted suspension over the entire area of a Mueller Hinton agar plate.
3. Apply the discs onto the inoculated plate, ensuring sufficient space between individual discs to allow for proper measurement of inhibition zones.
4. Incubate at $35 \pm 1^\circ\text{C}$ for 16-20 hours.

EVALUATING THE RESULTS

At the end of the incubation period, measure the inhibition halos and interpret as indicated in the following table.

Interpretative Table. CDT method for confirmation of ESBL in Enterobacteriaceae producing chromosomal AmpC.

Antibiotic Disc		Confirmation is positive if
Cefepime alone and with Clavulanic acid	FEP 30 μg and FEL 30+10 μg	≥ 5 mm increase in inhibition zone of cefepime with Clavulanic acid

QUALITY CONTROL

To check the performance of discs, MHA plates, inoculum and procedure used, test *Escherichia coli* ATCC® 25922 and *Escherichia coli* NCTC 13353^a (CTX-M-15 ESBL) as a negative and positive control for ESBL, respectively. The following criteria should be met:

Control strain	Criteria	Result
<i>Escherichia coli</i> NCTC 13353 ^a	≥ 5 mm increase in zone diameter with Cefepime (30 µg) + Clavulanic acid (10 µg) compared with Cefepime alone	Positive
<i>Escherichia coli</i> ATCC® 25922	± 1 mm difference in zone diameter with Cefepime (30 µg) + Clavulanic acid (10 µg) compared with Cefepime alone	Negative

^a CTX-M-15-producing strain *E. coli* NCTC 13353, which is recommended by CLSI for quality control of various cephalosporins and beta-lactam combination agents, is not included in the EUCAST guidelines; If discrete colonies or a haze of growth are present inside the zone of inhibition, measure the colony-free inner zone (CLSI M100 ED33:2023).

LIMITS

Diffusion susceptibility tests use an *in vitro* technique and cannot therefore reproduce the extremely complex *in vivo* conditions. Nevertheless, it is a useful and important tool that helps the clinician choose the correct therapy. Many variable factors influence the final result of the diffusion susceptibility test. The main ones are: the culture medium used, impregnation of the discs, inoculation of the medium, temperature, time and incubation atmosphere of the plates, pre-incubation and pre-diffusion conditions, depth of the medium, etc.

PRECAUTIONS

The disc cannot be classified as being hazardous according to current legislation but fall within the specific field of application where a safety data sheet must be supplied because they can cause phenomena of sensitization in sensitive subjects if they come into contact with the skin. The discs are disposable products. They are only for diagnostic *in vitro* use and are intended for professional use. They must be used in the laboratory by properly trained operators using approved aseptic and safety methods for pathogenic agents.

STORAGE

Store the unopened blister at -20°C to +8°C till the expiry date. Allow unopened cartridge to come to room temperature before removing it from the blister for minimising condensation on the discs. Leftover discs from an opened cartridge should be stored at 2-8°C for no more than 7 days. Return unused discs to the refrigerator as soon as the application of the discs has been completed. Dispose of expire discs.

ELIMINATING USED MATERIAL

After use, the discs and the material that comes into contact with the sample must be decontaminated and disposed of in accordance with current laboratory techniques for the decontamination and disposal of potentially infected material.

REFERENCES

- EUCAST guidelines for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance. Version 2.0, 2017.

Disc	µg	Disc	µg	Σ	REF
Cefepime FEP	30	Cefepime + Clavulanic acid FEL	40 (30+10)	100 Test	99006

TABLE OF SYMBOLS

LOT Batch code	IVD In Vitro Diagnostic Medical Device	Manufacturer	Use by	Fragile, handle with care
REF Catalogue number	Temperature limitation	Contains sufficient for <n> tests	Consult instructions for use	Do not reuse

Revision History

Revision	Release Date	Change Summary
1	2024-02-06	Revised TEST PROCEDURE (incubation) and QUALITY CONTROL; Updated REFERENCES

**LIOFILCHEM® s.r.l.**Via Scozia 64026 Roseto degli Abruzzi (TE) Italy
Tel. +39 0858930745 Fax +39 0858930330

www.liofilchem.com

liofilchem@liofilchem.com

ESBL (Chromos. Induc. AmpC) disc kit (acc. to EUCAST)

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