ANTIFUNGAL DISC

Discs for antifungal susceptibility test of yeasts

DESCRIPTION

ANTIFUNGAL DISC are paper discs with special features that are impregnated with antimicrobics, used for the susceptibility test of yeasts according to the CLSI method. Each antimicrobics is available in packages of 50 and 250 discs.

CONTENTS OF THE PACKAGES

The 50-test package contains 1 cartridge with 50 discs packaged in blister with silica gel.

The 250-test package contains 5 cartridges with 50 discs, packaged in blister with silica gel.

METHOD PRINCIPLE

The discs are applied to the surface of a culture medium inoculated with the suspension containing pure colonies of the microorganism under examination. After incubation, the plates are examined, the inhibition halos around each disc are examined and compared with the standard inhibition haloes: in this way the microorganisms are defined as being susceptible, susceptible-dose dependent, resistant or non susceptible to the tested antimicrobial agents.

COMPOSITION

The discs are made of high-quality paper in compliance with WHO and FDA specifications.

Each disc is impregnated with a standardised amount of antimicrobics.

GATHERING AND KEEPING SAMPLES

The colonies that are to be subjected to the susceptibility test are taken up by culture media that have been previously swabbed with the sample under examination. In the case of mixed colonies the strains must be purified in blood medium or Sabouraud Dextrose Agar before they are swabbed on the plates for the susceptibility test.

TEST PROCEDURE

For procedural details of the test, see the attached bibliography or the microbiology literature. A brief description of the CLSI method follows.

1. Take the blister of cartridges from the refrigerator and leave it for about an hour at ambient temperature in order to prevent condensation forming on the discs when it is opened as such condensation could affect long-term stability.
2. Swab 4 to 5 well isolated and morphologically similar colonies with a culture medium and suspend them in 5 ml of a sterile physiological solution until 0.5 McFarland turbidity is reached for *Candida* spp. and 1 McFarland for *C. neoformans*.
3. Immerse a sterile buffer in the suspension and squeeze it on the wall of the test tube to eliminate excess liquid. Drag it along the surface of the susceptibility test medium of yeasts (M.HINTON+GLUCOSE+METHYLEN BLUE AGAR, Ref. 11510).
4. Position the antifungal discs within 15 minutes from inoculation of the plates, pressing them with a sterile pliers on the surface of the agar and within 15 minutes of depositing of the discs leave the plate to incubate at 35±5°C for 24 hours; read at 48 hours only when insufficient growth is observed after 24 hours of incubation.

EVALUATING THE RESULTS

At the end of the incubation period, measure the inhibition halos and compare them with the diameters of the standard CLSI inhibition halos..

Pinpoint microcolonies at the zone edge or large colonies within a zone are encountered frequently and should be ignored.

CLINICAL INTERPRETATION

The susceptibility test carried out *in vitro* cannot exactly reproduce *in vivo* conditions. Nevertheless, it shows the effect of the concentration of the antimicrobics, which varies in the culture medium in relation to the growth of the microbial population.

The final choice of antimicrobics to administer to the patient is the responsibility of the clinician who possesses all the information on the patient.

QUALITY CONTROL

Each batch of **ANTIFUNGAL DISC** is subjected to precise and thorough checks in compliance with CLSI standards using the following strains: *Candida albicans* ATCC 90028; *Candida parapsilosis* ATCC 22019; *Candida tropicalis* ATCC 750; *Candida krusei* ATCC 6258.

LIMITS

Diffusion susceptibility tests use an *in vitro* technique and cannot therefore reproduce the extremely complex *in vivo* conditions. Nevertheless, it is a useful and important tool that helps the clinician choose the correct therapy. Many variable factors influence the final result of the diffusion susceptibility test. The main ones are: the culture medium used, impregnation of the discs, inoculation of the medium, temperature, time and incubation atmosphere of the plates, pre-incubation and pre-diffusion conditions, depth of the medium, etc.

PRECAUTIONS

The **ANTIFUNGAL DISC** cannot be classified as being hazardous according to current legislation. **ANTIFUNGAL DISC** are disposable products. **ANTIFUNGAL DISC** are for *in vitro* diagnostic use only and are intended for professional use. They must be used in the laboratory by properly trained operators using approved aseptic and safety methods for pathogenic agents.

STORAGE

Store the blister at -20 to 8°C till the expiry date. Leftover discs from an opened package should be stored at 2-8°C avoiding humidity until the expiry date. To further reduce humidity, store the discs inside the disc dispenser which contains desiccant. In such conditions, **ANTIFUNGAL DISC** can be used until the expiry date shown on the label. Do not use after this date. Dispose of if they show signs of deterioration.

ELIMINATING USED MATERIAL

After use, **ANTIFUNGAL DISC** and the material that comes into contact with the sample must be decontaminated and disposed of in accordance with current laboratory techniques for the decontamination and disposal of potentially infected material.

NOTE: See the **Antifungal Disc Interpretative Criteria and Quality Control** table at the end of this document.

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Antifungal Disc Interpretative Criteria and Quality Control.

ANTIFUNGAL AGENT	µg	Species	SDD				Zone Diameter Ranges (mm)		
			S ≥	I	SDD	R ≤	CLSI QUALITY CONTROL		
Caspofungin CAS	5	<i>C. albicans</i>	17	15-16	-	14	<i>C. albicans</i>	ATCC 90028	18-27
		<i>C. guilliermondii</i>	13	11-12	-	10	<i>C. parapsilosis</i>	ATCC 22019	14-23
		<i>C. krusei</i>	17	15-16	-	14	<i>C. tropicalis</i>	ATCC 750	20-27
		<i>C. parapsilosis</i>	13	11-12	-	10	<i>C. krusei</i>	ATCC 6258	19-26
		<i>C. tropicalis</i>	17	15-16	-	14			
Fluconazole FLU	25	<i>C. albicans</i>	17	-	14-16	13	<i>C. albicans</i>	ATCC 90028	28-39
		<i>C. glabrata</i>	-	-	≥ 15	14	<i>C. parapsilosis</i>	ATCC 22019	22-33
		<i>C. parapsilosis</i>	17	-	14-16	13	<i>C. tropicalis</i>	ATCC 750	26-37
		<i>C. tropicalis</i>	17	-	14-16	13			
Posaconazole POS	5		-	-	-	-	<i>C. albicans</i>	ATCC 90028	24-34
							<i>C. parapsilosis</i>	ATCC 22019	25-36
							<i>C. tropicalis</i>	ATCC 750	23-33
							<i>C. krusei</i>	ATCC 6258	23-31
Voriconazole VO	1	<i>C. albicans</i>	17	15-16	-	14	<i>C. albicans</i>	ATCC 90028	31-42
		<i>C. krusei</i>	15	13-14	-	12	<i>C. parapsilosis</i>	ATCC 22019	28-37
		<i>C. parapsilosis</i>	17	15-16	-	14	<i>C. krusei</i>	ATCC 6258	16-25
		<i>C. tropicalis</i>	17	15-16	-	14			

S = Susceptible I = Intermediate S-DD = Susceptible-Dose Dependent R = Resistant

PRESENTAZIONE/PRESENTATION/ERHÄLTICHE ANTIMYKOTIKA/PRESENTACIÓN

Description	CLSI	Ref.*	Description	CLSI	Ref.*
Amphotericin B AMB 20 µg		9071	Griseofulvin AGF 10 µg		9074
Amphotericin B AMB 10 µg		9137	Itraconazole ITC 50 µg		9107
Caspofungin CAS 5 µg ✓		9165 •	Itraconazole ITC 8 µg		9139
Clotrimazole CLO 50 µg		9097	Ketoconazole KCA 10 µg		9075
Econazole ECN 10 µg		9072	Ketoconazole KCA 15 µg		9140
Fluconazole FLU 100 µg		9069	Miconazole MCL 10 µg		9077
Fluconazole FLU 25 µg ✓		9166	Nystatin NY 100 units		9078
Flucytosine AFY 1 µg		9073	Posaconazole POS 5 µg ✓		9167
Flucytosine AFY 10 µg		9148	Voriconazole VO 1 µg ✓		9168

•, to be stored at -20°C

* Packaging

Each pack contains 5 cartridges of 50 discs (5x50 discs = 250 discs).

Single cartridges of 50 discs are available: add /1 to the catalogue ref. no. to indicate the relevant item.

Example: ref. 9168/1 indicates Voriconazole 1 µg in one single cartridge of 50 discs.

The above list of products might be out-of-date, [view the complete range of Antimicrobial Discs on Liofilchem's website](#)

BIBLIOGRAFIA/REFERENCES/BIBLIOGRAPHIE/REFERENCIAS

1. CLSI. Method for Antifungal Disk Diffusion Susceptibility Testing of Yeasts. 3rd ed. CLSI Guideline M44. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
2. CLSI. Performance Standards for Antifungal Susceptibility Testing of Yeasts. 1st ed. CLSI supplement M60. Wayne, PA: Clinical and Laboratory Standards Institute; 2017.

TABELLA DEI SIMBOLI/TABLE OF SYMBOLS/SYMBOLS/TABLA DE SÍMBOLOS

LOT Numero di Lotto Batch code Charge Código de lote	IVD Dispositivo Medico Diagnostico <i>in vitro</i> In vitro Diagnostic Medical Device Für die <i>in vitro</i> Diagnostik Dispositivo médico diagnóstico <i>in vitro</i>	 Fragile, maneggiare con cura Fragile, handle with care Zerbrechlich, vorsichtig behandeln Fragil, manipular con cuidado	 Non riutilizzare Do not reuse Einmalgebrauch, nicht wieder verwenden No reutilizar	 Fabricante Manufacturer Hersteller Fabricante
REF Numero di catalogo Catalogue number Artikelnummer Número de catálogo	 Limiti di Temperatura Temperature limitation Temperaturbereich für die Lagerung Límites de temperatura	 Contenuto sufficiente per <n> test Contains sufficient for <n> tests Enthält Material für <n> Tests Contenido suficiente para <n> análisis	 Attenzione vedere le istruzioni per l'uso Caution, consult accompanying documents Achtung, Packungsbeilage beachten Atención, consultar el documento adjunto	 Utilizzare entro Use by verwendbar bis Utilizar antes de

LIOFILCHEM® S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY
Tel +39 0858930745 Fax +39 0858930330 www.liofilchem.com

F00008