MULTODISCS

Anelli con otto antibiotici per test antibiogramma.

DESCRIZIONE

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MULTODISCS

Rings with eight antibiotics for susceptibility tests.

DESCRIPTION

MULTODISCS are rings of absorbent paper with eight arms, each terminating in a 6 mm disc impregnated with an antibiotic. Each disc is marked with an identification code. A barrier of inert material at the end of each arm prevents diffusion of the antibiotic. This makes it possible to have easily measurable circular inhibition zones and makes MULTODISCS suitable for sensitivity tests using Kirby-Bauer method (KB testing or disk diffusion antibiotic sensitivity testing). MULTODISCS are offered in a large variety of configurations.

CONTENTS OF THE PACKAGES

Each pack contains: 100 rings packed in a jar with a dryer pack and an instruction leaflet.

METHOD PRINCIPLE

The rings are applied to the surface of a culture medium inoculated with the culture broth prepared with pure colonies of the micro-organism under examination. After incubation the dishes are examined and the inhibition rings for each disc are measured and compared with the diameters of the standard inhibition ring. In this way the micro-organism can be shown as sensitive, intermediate or resistant to the antimicrobial agents tested.

COMPOSITION

The rings are made of high-quality paper in compliance with WHO and FDA specifications.

Antibiotic potency: the discs are manufactured according to DIN specification, i.e. the concentration of each antibiotic is within 90-125% of the concentration stated on the disc.

GATHERING AND KEEPING SAMPLES

The colonies that are to be subjected to the susceptibility testing are taken up by culture media that have been previously swabbed with the sample under examination. In the case of mixed colonies the bacterial strains must be purified before they are swabbed on the plates for the antibiogram.

TEST PROCEDURE

For procedural details of the test, see the attached bibliography or the microbiology literature. A brief description of the Kirby-Bauer method follows:

1. Take the jar from the refrigerator from the refrigerator and leave it for about an hour at ambient temperature in order to prevent condensation forming on the discs when it is opened as such condensation could affect long-term stability;
2. Swab 4 to 5 well isolated and morphologically similar colonies with a culture medium and suspend them in 5 ml of a suitable broth culture such as Tryptic Soy Broth;
3. Incubate the broth culture at 36 ± 1 °C until 0.5 McFarland turbidity is reached;
4. Immerse a sterile buffer in the broth culture or in a diluted form thereof and squeeze it on the wall of the test tube to eliminate excess liquid. Drag it along the surface of the antibiogram medium contained on the plate so as to produce even growth;
5. Position the rings within 15 minutes from inoculation of the plates, pressing them with a sterile pliers on the surface of the agar and within 15 minutes of depositing of the discs leave the plate to incubate according to the bacterial growth conditions.

EVALUATING THE RESULTS

At the end of the incubation period, measure the inhibition halos and interpret according to the current reference standards

CLINICAL INTERPRETATION

The susceptibility testing carried out *in vitro* cannot exactly reproduce *in vivo* conditions. Nevertheless, it shows the effect of the concentration of the antibiotic, which varies in the culture medium in relation to the growth of the microbial population.

The final choice of antibiotic to administer to the patient is the responsibility of the clinician who possesses all the information on the patient.

QUALITY CONTROL

Each batch of MULTODISCS is subjected to the quality control with the following bacterial strains:

<i>Enterococcus faecalis</i>	ATCC® 29212	<i>Staphylococcus aureus</i>	ATCC® 25923	<i>Haemophilus influenzae</i>	ATCC® 8468
<i>Escherichia coli</i>	ATCC® 25922	<i>Haemophilus influenzae</i>	ATCC® 49247	<i>Pseudomonas aeruginosa</i>	ATCC® 27853
<i>Staphylococcus aureus</i>	ATCC® 29213	<i>Haemophilus influenzae</i>	ATCC® 49766	<i>Streptococcus pneumoniae</i>	ATCC® 49619
<i>Candida albicans</i> *	ATCC® 90028	<i>Candida krusei</i> *	ATCC® 6258	<i>Candida parapsilosis</i> *	ATCC® 22019
<i>Candida tropicalis</i> *	ATCC® 750	* only for "Yeasts"			

LIMITS

Diffusion antibiograms use an *in vitro* technique and cannot therefore reproduce the extremely complex *in vivo* conditions. Nevertheless, it is a useful and important tool that helps the clinician choose the correct therapy. Many variable factors influence the final result of the diffusion antibiogram. The main ones are: the culture medium used, impregnation of the discs, inoculation of the medium, temperature, time and incubation atmosphere of the plates, pre-incubation and pre-diffusion conditions, depth of the medium, etc.

PRECAUTIONS

The MULTODISCS cannot be classified as being hazardous according to current legislation but fall within the specific field of application where a safety data sheet must be supplied because they can cause phenomena of sensitization in sensitive subjects if they come into contact with the skin. MULTODISCS are disposable products. MULTODISCS are only for diagnostic *in vitro* use and are intended for professional use. They must be used in the laboratory by properly trained operators using approved aseptic and safety methods for pathogenic agents.

STORAGE

MULTODISCS must be stored for long periods of time at -20 °C in their original packaging. For short periods, MULTODISCS can be stored at 2-8 °C, provided that they are still in their original packaging. Do not store them near sources of heat and do not expose them to excessive temperature variations. The packages must be placed in the refrigerator as soon as possible after use. In such conditions, MULTODISCS can be used until the expiry date shown on the label. Do not use after this date. Dispose of if they show signs of deterioration.

ELIMINATING USED MATERIAL

After use, MULTODISCS and the material that comes into contact with the sample must be decontaminated and disposed of in accordance with current laboratory techniques for the decontamination and disposal of potentially infected material.

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1. CLSI M100-S23. Performance Standards for Antimicrobial Susceptibility Testing, January 2013.
2. CLSI M02-A10. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard-Tenth Edition.
3. CLSI document M44-A2—Method for Antifungal Disk Diffusion Susceptibility Testing of Yeasts; Approved Guideline—Second Edition (August 2009).
4. CLSI document M44-S3—Zone Diameter Interpretive Standards, Corresponding Minimal Inhibitory Concentration (MIC) Interpretive Breakpoints, and Quality Control Limits for Antifungal Disk Diffusion Susceptibility Testing of Yeasts; Third Informational Supplement.
5. EUCAST. Breakpoint tables for interpretation of MICs and zone diameters. Version 3.0, January 2013.
6. FDA (1978) Codes of Fed.Rebs. 21.Part 460.
7. WHO (1977) Tech Rep.Ser.no.610.
8. DIN 58940-2 Medical microbiology - Susceptibility testing of microbial pathogens to antimicrobial agents - Part 2: Active substance carriers for the agar diffusion test; 2007-10.


**TABELLA DEI SIMBOLI / TABLE OF SYMBOLS / TABLA DE LOS SÍMBOLOS
TABLEAU DES SYMBOLE / ΠΙΝΑΚΑΣ ΣΥΜΒΟΛΩ / MESA DE SÍMBOLOS**

IVD	Dispositivo medico-diagnostico <i>in vitro</i> / <i>In Vitro</i> Diagnostic Medical Device Producto sanitario para diagnóstico <i>in vitro</i> / Dispositif médical de diagnostic <i>in vitro</i> <i>In Vitro</i> Διαγνωστικό Ιατροτεχνολογικό προϊόν / Dispositivo médico para diagnóstico <i>in vitro</i>
	Non riutilizzare / Do not reuse / No reutilizar Ne pas réutiliser / Μην κάνετε επαναληπτική χρήση / Não reutilizar
	Fabbricante / Manufacturer / Fabricante Fabricant / Κατασκευαστής / Fabricante
	Contenuto sufficiente per "n" saggi / Contains sufficient for <n> tests / Contenido suficiente para <n> ensayos Contenu suffisant pour "n" tests / Περιεχόμενο επαρκές για «n» εξετάσεις / Conteúdo suficiente para "n" ensaios
REF	Numero di catalogo / Catalogue number / Número de catálogo Référence du catalogue / Αριθμός καταλόγου / Referência de catálogo
	Fragile, maneggiare con cura / Fragile, handle with care / Frágil, manipular con precaución Fragile, manipuler avec précautions / Εύθραυστο, να χρησιμοποιείται με προσοχή / Frágil, manusear com cuidado
	Utilizzare entro / Use By / Fecha de caducidad Utiliser jusque / Ημερομηνία λήξης / Prazo de validade
	Attenzione, vedere le istruzioni per l'uso / Caution, consult accompanying documents Atención, ver instrucciones de uso / Attention voir notice d'instructions Προειδοποίηση, συμβουλευτείτε τα συνοδά έντυπα / Atenção, consulte a documentação incluída
	Limiti di temperatura / Temperature limitation / Limite de temperatura Limites de température / Περιορισμοί θερμοκρασίας / Limites de temperatura
LOT	Codice del lotto / Batch code / Código de lote Code du lot / Αριθμός Παρτίδας / Código do lote

**LIOFILCHEM® s.r.l.**

Via Scozia Zona Ind.le - 64026 Roseto D.A. (TE) - Italy

Tel. +390858930745 Fax +390858930330 Website: www.liofilchem.net E-Mail: liofilchem@liofilchem.net

F01731
Rev.1 / 02.04.2015