SirScan Discs

Antibiotic discs for susceptibility tests

Contents	Page
English	1
Italiano	2
Español	3
Deutsch	4
Français	5
Product list	7
References	7

SirScan Discs

Antibiotic discs for susceptibility tests

ENGLISH

DESCRIPTION

SirScan Discs are paper discs with special features, that are impregnated with antibiotic and used for the susceptibility test according to the Kirby-Bauer antibiotic testing (KB testing or disk diffusion antibiotic sensitivity testing).

CONTENTS OF THE PACKAGES

The 250-test box contains 5 cartridges of 50 discs, each cartridge individually packed in a desiccant envelope.

METHOD PRINCIPLE

The discs are applied to the surface of a culture medium inoculated with a pure colony suspension of the microorganism under examination. After incubation, the plates are examined, the inhibition halos around each disc are examined and compared with the standard inhibition haloes: in this way the microorganisms are defined as being susceptible, intermediate or resistant to the tested antimicrobial agents.

COMPOSITION

SirScan Discs for antimicrobial susceptibility test are made of high-quality paper specifications.

The discs are manufactured under the quality systems UNI EN ISO 9001:2015 and EN ISO 13485:2016. The concentration of each antibiotic is within 90-125% of the concentration stated on the disc.

GATHERING AND KEEPING SAMPLES

The colonies that are to be subjected to the susceptibility test are taken up by culture media that have been previously swabbed with the sample under examination. In the case of mixed colonies the bacterial strains must be purified before they are swabbed on the plates for the susceptibility test.

TEST PROCEDURE

1. Allow discs to equilibrate to room temperature before opening the container (cartridge or canister) in order to minimize the condensation on the discs, which could affect long-term stability.
2. Make a suspension of the test organism following the chosen method (0.5 McF or Inoclic).
3. Using a swab, inoculate a suitable agar plate medium by uniformly spreading the suspension over the entire agar surface.
4. Apply discs firmly to the surface of the inoculated agar plate.
5. Incubate plates at the appropriate temperature and conditions following the reference methodology (CLSI, EUCAST).
6. As soon as the application of the discs has been completed, pack away unused discs following supplier's storage indications.

NOTE 1: The medium to be used depends on the organism under investigation and the methodology followed, and must be validated by the media manufacturer for antimicrobial susceptibility testing.

NOTE 2: It is recommended to apply the 3x15 minutes rule, to carry out the susceptibility test:

-Use the inoculum suspension within 15 minutes of preparation;

-Apply discs within 15 minutes of inoculation;

-Incubate plates within 15 minutes of disc application.

For more details, please refer to the current published standards.

EVALUATING THE RESULTS / QUALITY CONTROL

At the end of the incubation period, measure the inhibition halos and interpret according to the current reference standards (CLSI, EUCAST).

CLINICAL INTERPRETATION

The susceptibility test carried out *in vitro* cannot exactly reproduce *in vivo* conditions. Nevertheless, it shows the effect of the concentration of the antibiotic, which varies in the culture medium in relation to the growth of the microbial population.

The final choice of antibiotic to administer to the patient is the responsibility of the clinician who possesses all the information on the patient.

LIMITS

Diffusion susceptibility tests use an *in vitro* technique and cannot therefore reproduce the extremely complex *in vivo* conditions. Nevertheless, it is a useful and important tool that helps the clinician choose the correct therapy. Many variable factors influence the final result of the diffusion susceptibility test. The main ones are: the culture medium used, impregnation of the discs, inoculation of the medium, temperature, time and incubation atmosphere of the plates, pre-incubation and pre-diffusion conditions, depth of the medium, etc.

PRECAUTIONS

SirScan Discs cannot be classified as being hazardous according to current legislation. SirScan Discs are disposable products. SirScan Discs are only for diagnostic *in vitro* use and are intended for professional use. They must be used in the laboratory by properly trained operators using approved aseptic and safety methods for pathogenic agents.

STORAGE

The unopened package of Antibiotic Disc can be stored in most cases at -20°C to $+8^{\circ}\text{C}$ till the expiry date. Some products have to be stored at -20°C as maximum storage temperature. The recommended temperature limits can be found both on the product envelope and on the box label. Leftover discs from an opened CARTRIDGE need to be stored at $2-8^{\circ}\text{C}$ for no more than 7 days.

ELIMINATING USED MATERIAL

After use, SirScan Discs and the material that comes into contact with the sample must be decontaminated and disposed of in accordance with current laboratory techniques for the decontamination and disposal of potentially infected material.

VIGILANCE

Any serious incident occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user is established.

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•, to be stored at -20°C

Storage Unless otherwise indicated, discs may be stored at any temperature between -20°C and +8°C.

CE Marking Unless otherwise indicated, all products intended for clinical applications are CE marked according to the European Directive 98/79/EC for In Vitro Diagnostic Medical Devices.

REFERENCES

1. The European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MICs and zone diameters. Version 11.0, 2021. <http://www.eucast.org>
2. The European Committee on Antimicrobial Susceptibility Testing. Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST, Version 11.0, 2021. <http://www.eucast.org>
3. EUCAST Disk Diffusion Method for Antimicrobial Susceptibility Testing - Version 9.0 (January 2021).
4. EUCAST technical guidance on the use of the combination disk test (CDT) for confirmation of ESBL in Enterobacterales. New disk potencies for combination disks containing cefotaxime and ceftazidime without and with clavulanic acid (February 2019).
5. EUCAST guidelines for detection of resistance mechanisms and specific resistances of clinical and/or epidemiological importance. Version 2.0 (July 2017).
6. DIN 58940-2 Medical microbiology - Susceptibility testing of microbial pathogens to antimicrobial agents - Part 2: Active substance carriers for the agar diffusion test; 2007-10.
7. FDA (1978) Codes of Fed.Rebs. 21.Part 460.
8. WHO (1977) Tech rep.Ser.n°610.

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