


Disc kits for Antimicrobial Resistance Testing

KPC&MBL&OXA-48 discs kit


Description		µg	Packaging	Ref.
Meropenem	MRP	10	50 Test	99007
Meropenem + Phenylboronic acid	MR+BO			
Meropenem + Cloxacillin	MR+CL			
Meropenem + EDTA	MR+EDTA			
Temocillin	TMO	30		



Klebsiella pneumoniae
clinical isolate
OXA-48 positive

KPC&MBL discs kit (according to EUCAST)

Description		µg	Packaging	Ref.
Meropenem	MRP	10	50 Test	99003
Meropenem + Phenylboronic acid	MR+BO			
Meropenem + Cloxacillin	MR+CL			
Meropenem + EDTA	MR+EDTA			



Klebsiella pneumoniae
ATCC® BAA-1705™
KPC positive

EUCAST ESBL discs kit

Description		µg	Packaging	Ref.
Cefotaxime	CTX	5	50 Test	99002
Cefotaxime + Clavulanic acid	CTL	15 (5+10)		
Ceftazidime	CAZ	10		
Ceftazidime + Clavulanic acid	CAL	20 (10+10)		
Cefepime	FEP	30		
Cefepime + Clavulanic acid	FEL	40 (30+10)		



Liofilchem srl

Zona Industriale, 64026 Roseto degli Abruzzi (TE), Italy

Tel +39 0858930745

www.liofilchem.com

Headquarter, Manufacturing Site, International distribution

In US available as RUO device

Liofilchem, Inc.

465 Waverley Oaks Rd. Suite 317

Waltham, MA 02452, USA


Tel 781-902-0312

US Distribution Center

Disc kits for Antimicrobial Resistance Testing

ESBL discs kit (according to CLSI)


Description		µg	Packaging	Ref.
Cefotaxime	CTX	30	50 Test	99005
Cefotaxime + Clavulanic acid	CTL	40 (30+10)		
Ceftazidime	CAZ	30		
Ceftazidime + Clavulanic acid	CAL	40 (30+10)		



Klebsiella pneumoniae subsp. *pneumoniae*
ATCC® 700603™
ESBL positive

ESBL (chromosomal ind. AmpC) discs kit (according to EUCAST)

Description		µg	Packaging	Ref.
Cefepime	FEP	30	100 Test	99006
Cefepime + Clavulanic acid	FEL	40 (30+10)		



Enterobacter cloacae
clinical isolate
ESBL positive

ESBL+AmpC screen discs kit

Description	µg	Packaging	Ref.
Cefotaxime	30	50 Test	99008
Cefotaxime + Clavulanic acid	40 (30+10)		
Cefotaxime + Cloxacillin			
Cefotaxime + Clavulanic acid + Cloxacillin			

AmpC discs kit

Description	µg	Packaging	Ref.
Cefotaxime	30	50 Test	99009
Cefotaxime + Cloxacillin			
Ceftazidime	30		
Ceftazidime + Cloxacillin			



Liofilchem srl

Zona Industriale, 64026 Roseto degli Abruzzi (TE), Italy

Tel +39 0858930745

www.liofilchem.com

Headquarter, Manufacturing Site, International distribution

In US available as RUO device

Liofilchem, Inc.

465 Waverley Oaks Rd. Suite 317

Waltham, MA 02452, USA

Tel 781-902-0312

US Distribution Center