



# Rotating Plate Inoculator



for confluent microbial growth



The Rotating Plate Inoculator (RPI) is used to inoculate 90 mm agar plates for AST methods such as MIC Test Strip and antibiotic disks.

The rotation speed can be easily adjusted to user preference.

The RPI allows for standardization and even distribution of the inoculum.



**LIOFILCHEM® s.r.l.**

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY







Tel +39 0858930745 Fax +39 0858930330 Website: [www.liofilchem.net](http://www.liofilchem.net) E-mail: [liofilchem@liofilchem.net](mailto:liofilchem@liofilchem.net)

# Rotating Plate Inoculator



Ref. 96794

**Ideal for manual agar-based  
antimicrobial susceptibility testing methods**



**LIOFILCHEM® s.r.l.**

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY

Tel +39 0858930745 Fax +39 0858930330 Website: [www.liofilchem.net](http://www.liofilchem.net) E-mail: [liofilchem@liofilchem.net](mailto:liofilchem@liofilchem.net)

Rev.0 / 08.09.2017