

Step 1

- Input patient and specimen information (LIS/LIMS or manual)
- Insert test plate
- Select "Read AST Plate" for instant color digital reading

Productivity Enhancement Features:
 Bi-Directional LIS/LIMS
 21.5 inch Touch Screen Monitor
 Bar Code Reader
 Rapid-Batch Entry Mode




Step 2

- Review plate image and zone readings
- Select "View Results"

User may adjust zone diameter on-screen.
 Plate images can be saved and archived.

Touch Screen



Step 3

- Review test results
- Print and transfer test results to LIS/LIMS



Epidemiology
 Reports and
 Charts

Print and transfer test results to LIS/LIMS

