

## Insert Identification Panel Select "Read ID Panel"




- API<sup>®</sup> (bioMérieux)
- RapID<sup>™</sup> (Remel)
- Crystal<sup>™</sup> (BD)

Read and interpret all varieties of ID panels




21.5 inch  
Touch Screen

## View Results



The BIOMIC V3 provides:  
Instant reading of +/- reactions  
Automatic organism identification look-up  
Saved ID panel color images


## Additional Identification Functions



96-Well Microtiter Reading  
(optional feature)



Chromogenic Agar Reading



IDTree<sup>™</sup> Software  
Provides training based on  
M35 CLSI methods