


## Instant Reading



QC AST plates and ID panels are instantly read and interpreted. All images can be saved and archived.



## QC Test Results

Drug Name	ZD (mm)	ZD Limits (mm)	Quality
Ampicillin	30	27 - 35	OK
Clindamycin	27	24 - 30	OK
Cefotaxime	26	25 - 31	OK
Cefazolin	29	29 - 35	OK
Erythromycin	28	22 - 30	OK
Imipenem	39		N/A
Levofloxacin	30	25 - 30	OK
Oxacillin	21	18 - 24	OK
Penicillin G	31	28 - 37	OK
Trimethoprim-sulfam	29	24 - 32	OK
Tetracycline	29	24 - 30	OK
Vancomycin	18	17 - 21	OK

QC results are automatically interpreted and saved following CLSI, EUCAST or BSAC guidelines.

## QC Summary Reports



QC summary reports are generated in a variety of formats. This provides electronic documentation and facilitates laboratory inspections.