

## ComASP® Cefiderocol 0.008-128 µg/mL

Broth microdilution panel for Cefiderocol susceptibility testing according to ISO 20776-1



ComASP® Cefiderocol 0.008-128 µg/mL is a 2-test panel containing the dried-up antibiotic in 15 two-fold dilutions for the quantitative determination of the minimum inhibitory concentration (MIC) of Cefiderocol against Gram-negative non-fastidious organisms such *Enterobacterales* and *Pseudomonas aeruginosa*, *Acinetobacter* spp., *Stenotrophomonas maltophilia*.

Kit content:

- 4 individually packed ComASP® Cefiderocol 0.008-128 µg/mL panels (2 samples/panel)
- 8 Tubes of **Iron-Depleted Cation-adjusted Mueller Hinton Broth (ID-CAMHB)** (3.6 ml)
- Sealing Film

Ref.	Device format	Description	Packaging
75009	ComASP®	Cefiderocol 0.008-128 µg/mL	8 tests



**Liofilchem srl**

Zona Industriale, 64026 Roseto degli Abruzzi (TE), Italy

Tel +39 0858930745

[www.liofilchem.com](http://www.liofilchem.com)

Headquarters, Manufacturing Site, International distribution


**Liofilchem, Inc.**

465 Waverley Oaks Rd. Suite 317

Waltham, MA 02452, USA

Tel 781-902-0312

US Distribution Center



**MTS™ Cefiderocol**  
**0.016-256 µg/mL**  
*for P. aeruginosa*



**Cefiderocol disc**  
**30 µg**

Ref.	Device format	Description	Packaging
92067	MTS™ (MIC Test Strip)	Cefiderocol 0.016-256 µg/mL	30 mic test
920670	MTS™ (MIC Test Strip)	Cefiderocol 0.016-256 µg/mL	100 mic test
920671	MTS™ (MIC Test Strip)	Cefiderocol 0.016-256 µg/mL	10 mic test
9266	Antibiotic disc in cartridge	Cefiderocol 30 µg	5x50 Discs
9266/1	Antibiotic disc in cartridge	Cefiderocol 30 µg	50 Discs
9266/2	Antibiotic disc in canister	Cefiderocol 30 µg	250 Discs



**Liofilchem srl**

Zona Industriale, 64026 Roseto degli Abruzzi (TE), Italy  
 Tel +39 0858930745  
 www.liofilchem.com  
 Headquarters, Manufacturing Site, International distribution

**Liofilchem, Inc.**

465 Waverley Oaks Rd. Suite 317  
 Waltham, MA 02452, USA  
 Tel 781-902-0312  
 US Distribution Center