


Liofilchem® CONTACT SLIDE

flexible contact slide for microbiological screening of surfaces

ISO 18593



Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.

printed in Italy

© Liofilchem 2020

Rev. 7 - 28.05.2020

ref. 6550616

Description	Packaging	Ref.	# Product with minimum order
CONTACT SLIDE Chrom 1 (Chromatic™ E.coli-coliform / Chromatic™ Salmonella) Chromogenic selective media for detecting <i>E. coli</i> , coliforms and <i>Salmonella</i> spp. on surfaces or in food and waters.	20 slides 120 slides	525292 # 53529 #	
CONTACT SLIDE Chrom 2 (Chromatic™ E.coli-coliform / PCA+TTC + Neutralizing) Chromogenic selective medium for detecting <i>E. coli</i> , coliforms and for total bacterial count.	20 slides	525392	
CONTACT SLIDE Chrom 3 (Chromatic™ E.coli-coliform / Baird Parker Agar + Neutralizing) Chromogenic selective medium for detecting <i>E. coli</i> , coliforms and <i>Staphylococcus aureus</i> with inactivation of disinfectants.	20 slides	525462 #	
CONTACT SLIDE 1 (PCA + TTC + Neutralizing / V.R.B.G. Agar + Neutralizing) Total bacterial count, enterobacteria, with inactivation of disinfectants.	20 slides 120 slides	525262 53526	
CONTACT SLIDE 2 (PCA + TTC + Neutralizing / Rose Bengal CAF Agar + Neutralizing) Total bacterial count, yeasts and moulds, with inactivation of disinfectants.	20 slides 120 slides	525272 53527	
CONTACT SLIDE 3 (PCA+TTC + Neutralizing / Vogel-Johnson Agar) Total bacterial count with inactivation of disinfectants, staphylococci.	20 slides 120 slides	525302 53530	
CONTACT SLIDE 4 (Cetrimide Agar + Neutralizing / Rose Bengal CAF Agar + Neutralizing) <i>Pseudomonas</i> spp., yeasts and moulds.	20 slides 120 slides	525282 # 53528 #	
CONTACT SLIDE 5 (V.R.B.G. Agar + Neutralizing / Bile Esculin Agar + Neutralizing) Enterobacteria, fecal streptococci, with inactivation of disinfectants.	20 slides 120 slides	525312 53531	
CONTACT SLIDE 8 (PCA + TTC + Neutralizing / V.R.B.L. Agar + Neutralizing) Total bacterial count, coliforms, with inactivation of disinfectants.	20 slides 120 slides	525622 53562	
CONTACT SLIDE 9 (PCA + TTC + Neutralizing / PCA + TTC + Neutralizing) Total bacterial count with inactivation of disinfectants.	20 slides 120 slides	525322 53532	
CONTACT SLIDE 10 (Rose Bengal CAF Agar + Neutralizing / Vogel-Johnson Agar) Yeasts and moulds with inactivation of disinfectants, staphylococci.	20 slides 120 slides	525452 # 53545 #	
CONTACT SLIDE TSA + Neutralizing / TSA + Neutralizing Total bacterial count with inactivation of disinfectants.	20 slides	525472	
CONTACT SLIDE TSA + Neutralizing / Rose Bengal CAF Agar + Neutralizing Total bacterial count with inactivation of disinfectants, yeasts and moulds.	20 slides 120 slides	525482 53548	
CONTACT SLIDE PCA + TTC + Neutralizing / Dichloran Rose Bengal CAF + Neutralizing Total bacterial count with inactivation of disinfectants, yeasts and moulds.	20 slides 120 slides	525432 # 53543 #	



Liofilchem srl

Zona Industriale, 64026 Roseto degli Abruzzi (TE), Italy
Tel +39 0858930745
www.liofilchem.com

Headquarters, Manufacturing Site, International distribution

Liofilchem, Inc.

465 Waverley Oaks Rd. Suite 317
Waltham, MA 02452, USA
Tel 781-902-0312

US Distribution Center



Microbiology products since 1983