


Liofilchem®

# CONTACT SLIDE

flexible contact slide for microbiological screening of surfaces

ISO 18593



Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.

Description	Packaging	Ref.	# Product with minimum order
<b>CONTACT SLIDE Chrom 1</b> (Chromatic™ E.coli-coliform / Chromatic™ Salmonella) Chromogenic selective media for detecting <i>E. coli</i> , coliforms and <i>Salmonella</i> spp. on surfaces or in food and waters.	20 slides 120 slides	525292 # 53529 #	
<b>CONTACT SLIDE Chrom 2</b> (Chromatic™ E.coli-coliform / PCA+TTC + NEUTRALIZING) Chromogenic selective medium for detecting <i>E. coli</i> , coliforms and for total bacterial count.	20 slides 120 slides	525392 53539	
<b>CONTACT SLIDE Chrom 3</b> (Chromatic™ E.coli-coliform / BAIRD PARKER AGAR + NEUTRALIZING) Chromogenic selective medium for detecting <i>E. coli</i> , coliforms and <i>Staphylococcus aureus</i> with inactivation of disinfectants.	20 slides 120 slides	525462 # 53546 #	
<b>CONTACT SLIDE 1</b> (PCA + TTC + NEUTRALIZING / V.R.B.G. AGAR + NEUTRALIZING) Total bacterial count, enterobacteria, with inactivation of disinfectants.	20 slides 120 slides	525262 53526	
<b>CONTACT SLIDE 2</b> (PCA + TTC + NEUTRALIZING / ROSE BENGAL CAF AGAR + NEUTRALIZING) Total bacterial count, yeasts and moulds, with inactivation of disinfectants.	20 slides 120 slides	525272 53527	
<b>CONTACT SLIDE 3</b> (PCA+TTC + NEUTRALIZING / VOGEL JOHNSON AGAR) Total bacterial count with inactivation of disinfectants, staphylococci.	20 slides 120 slides	525302 53530	
<b>CONTACT SLIDE 4</b> (CETRIMIDE AGAR + NEUTRALIZING / ROSE BENGAL CAF AGAR + NEUTRALIZING) <i>Pseudomonas</i> spp., yeasts and moulds.	20 slides 120 slides	525282 # 53528 #	
<b>CONTACT SLIDE 5</b> (V.R.B.G. AGAR + NEUTRALIZING / BILE ESCULIN AGAR + NEUTRALIZING) Enterobacteria, fecal streptococci, with inactivation of disinfectants.	20 slides 120 slides	525312 53531	
<b>CONTACT SLIDE 8</b> (PCA + TTC + NEUTRALIZING / V.R.B.L. AGAR + NEUTRALIZING) Total bacterial count, coliforms, with inactivation of disinfectants.	20 slides 120 slides	525622 53562	
<b>CONTACT SLIDE 9</b> (PCA + TTC + NEUTRALIZING / PCA + TTC + NEUTRALIZING) Total bacterial count with inactivation of disinfectants.	20 slides 120 slides	525322 53532	
<b>CONTACT SLIDE 10</b> (ROSE BENGAL CAF AGAR + NEUTRALIZING / VOGEL JOHNSON AGAR) Yeasts and moulds with inactivation of disinfectants, staphylococci.	20 slides 120 slides	525452 # 53545 #	
<b>CONTACT SLIDE TSA + NEUTRALIZING</b> / TSA + NEUTRALIZING Total bacterial count with inactivation of disinfectants.	20 slides 120 slides	525472 # 53547 #	
<b>CONTACT SLIDE TSA + NEUTRALIZING</b> / ROSE BENGAL CAF AGAR + NEUTRALIZING Total bacterial count with inactivation of disinfectants, yeasts and moulds.	20 slides 120 slides	525482 # 53548 #	



**Liofilchem s.r.l.**  
Zona Industriale  
64026, Roseto degli Abruzzi (Te), Italy  
tel +39 0858930745 fax +39 0858930330

[www.liofilchem.net](http://www.liofilchem.net)

ref. 6550616  
Rev. 5.1 - 30.06.2014  
© Liofilchem 2016  
printed in Italy

# Liofilchem® DIP-SLIDE

## Dip-Slide for the microbiological screening in **Urinary Tract Infections**



Description	Packaging	Ref.
<b>URITEST</b> (MAC CONKEY AGAR / CETRIMIDE AGAR / CLED ANDRADE AGAR) Enterobacteria, <i>Pseudomonas</i> spp., total bacterial count.	20 slides 120 slides	500152 51015
<b>URITEST N</b> (MAC CONKEY AGAR / CETRIMIDE AGAR / CLED AGAR) Enterobacteria, <i>Pseudomonas</i> spp., total bacterial count.	20 slides 120 slides 500 slides	500232 51023 51123
<b>URITEST C</b> (MAC CONKEY MUG AGAR / CETRIMIDE AGAR / CLED ANDRADE AGAR) Enterobacteria, <i>Escherichia coli</i> , <i>Pseudomonas</i> spp., total bacterial count.	20 slides 120 slides	500242 51024
<b>URITEST EC</b> (MAC CONKEY AGAR / E.COLI AGAR / CLED AGAR) Enterobacteria, <i>Escherichia coli</i> and coliforms, total bacterial count.	20 slides 120 slides	500412 51041
<b>URITEST 2</b> (CLED AGAR/ MAC CONKEY AGAR) Total bacterial count, enterobacteria.	20 slides 120 slides 500 slides	500302 51030 51130
<b>URITEST M</b> (MAC CONKEY AGAR / MALT SELECTIVE AGAR / CLED ANDRADE AGAR) Enterobacteria, yeasts, total bacterial count.	20 slides 120 slides 500 slides	500182 51018 51118
<b>URITEST MALTO</b> (MAC CONKEY AGAR / MALT AGAR / CLED AGAR) Enterobacteria, yeasts, total bacterial count.	20 slides 120 slides 500 slides	500402 51040 51140
<b>URITEST EF</b> (MAC CONKEY AGAR / BILE ESCULIN AGAR / CLED AGAR) Enterobacteria, enterococci, total bacterial count.	20 slides 120 slides	500702 51070
<b>URITEST PENTA</b> (CLED AGAR / MAC CONKEY AGAR / CETRIMIDE AGAR / MALT AGAR / BILE ESCULIN AGAR) Total bacterial count, enterobacteria, <i>Pseudomonas</i> spp., yeasts, enterococci.	20 slides 120 slides	500142 51014

## Dip-Slide for the microbiological screening in **Urogenital Infections**

Description	Packaging	Ref.
<b>VAGITEST</b> (CHOCOLATE / ROGOSA / THAYER MARTIN) Fastidious microorganisms, lactobacilli, <i>Neisseria gonorrhoeae</i> .	10 slides 120 slides	50020 51020

## Dip-Slide for the microbiological screening in **Skin Infections**

Description	Packaging	Ref.
<b>DERMATEST</b> (MALTO / CORN MEAL / DTM-TAPLIN) Pathogenic dermatophytes and others non pathogenic fungi.	10 slides 120 slides	50021 51021

Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.

[www.liofilchem.net](http://www.liofilchem.net)



**Liofilchem s.r.l.**  
Zona Industriale  
64026, Roseto degli Abruzzi (Te), Italy  
tel +39 0858930745 fax +39 0858930330

ref. 6550616  
Rev. 5.1 - 30.06.2014  
© Liofilchem 2016  
printed in Italy