

Liofilchem® DIP-SLIDE

Dip-Slide for the microbiological screening in Urinary Tract Infections



Description	Packaging	Ref.	# Product with minimum order
URITEST (MacConkey Agar / Cetrimide Agar / CLED Andrade Agar) Enterobacteria, <i>Pseudomonas</i> spp., total bacterial count.	20 slides	500152	
	120 slides	51015	
URITEST N (MacConkey Agar / Cetrimide Agar / CLED Agar) Enterobacteria, <i>Pseudomonas</i> spp., total bacterial count.	20 slides	500232	
	120 slides	51023	
	500 slides	51123	
URITEST C (MacConkey MUG Agar / Cetrimide Agar / CLED Andrade Agar) Enterobacteria, <i>Escherichia coli</i> , <i>Pseudomonas</i> spp., total bacterial count.	20 slides	500242 #	
	120 slides	51024 #	
URITEST EC (MacConkey Agar / E.coli Agar / CLED Agar) Enterobacteria, <i>Escherichia coli</i> and coliforms, total bacterial count.	20 slides	500412	
	120 slides	51041	
URITEST 2 (CLED Agar/ MacConkey Agar) Total bacterial count, enterobacteria.	20 slides	500302	
	120 slides	51030	
URITEST M (MacConkey Agar / Malt Selective Agar / CLED Andrade Agar) Enterobacteria, yeasts, total bacterial count.	20 slides	500182	
	120 slides	51018	
URITEST MALTO (MacConkey Agar / Malt Agar / CLED Agar) Enterobacteria, yeasts, total bacterial count.	20 slides	500402	
	120 slides	51040	
URITEST EF (MacConkey Agar / Bile Esculin Agar / CLED Agar) Enterobacteria, enterococci, total bacterial count.	20 slides	500702	
	120 slides	51070	
URITEST PENTA (CLED Agar / MacConkey Agar / Cetrimide Agar / Malt Agar / Bile Esculin Agar) Total bacterial count, enterobacteria, <i>Pseudomonas</i> spp., yeasts, enterococci.	20 slides	500142	
	120 slides	51014	

Dip-Slide for the microbiological screening in Urogenital Infections

Description	Packaging	Ref.
VAGITEST (Chocolate / Rogosa / Thayer Martin) Fastidious microorganisms, lactobacilli, <i>Neisseria gonorrhoeae</i> .	10 slides	50020

Dip-Slide for the microbiological screening in Skin Infections

Description	Packaging	Ref.
DERMATEST (Malto / Corn Meal / DTM-TAPLIN) Pathogenic dermatophytes and others non pathogenic fungi.	10 slides	50021

Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.

printed in Italy

© Liofilchem 2021

Rev. 8 - 01.03.2021

ref. 6550616



CE IVD

Liofilchem srl

Zona Industriale, 64026 Roseto degli Abruzzi (TE), Italy

Tel +39 0858930745

www.liofilchem.com

Headquarters, Manufacturing Site, International distribution

In US available as RUO device

Liofilchem, Inc.

465 Waverley Oaks Rd. Suite 317

Waltham, MA 02452, USA

Tel 781-902-0312

US Distribution Center



Microbiology products since 1983