


Multodisc with eight antimicrobial agents

*Customized
 antibiotics &
 concentrations
 welcome!*



page 1 of 2


Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.

Description		µg	CLSI	EUCAST	Packaging	Ref.
ANAEROBES	Carbenicillin	CAR 100			100 Multodiscs	95200
	Cefoxitin	FOX 30				
	Clindamycin	CD 2				
	Chloramphenicol	C 30				
	Erythromycin	E 15				
	Metronidazole	MTZ 5				
	Ampicillin	AMP 10				
	Tetracycline	TE 30				
ENTEROCOCCI	Tigecycline	TGC 15		✓	100 Multodiscs	95210
	Ampicillin	AMP 2		✓		
	Linezolid	LNZ 10		✓		
	Nitrofurantoin	F 100		✓		
	Trimethoprim-Sulfamethoxazole	SXT 25	✓	✓		
	Vancomycin	VA 5		✓		
	Gentamicin	CN 30		✓		
	Teicoplanin	TEC 30	✓	✓		
ENTEROBACTERIA URINE	Levofloxacin	LEV 5	✓	✓	100 Multodiscs	95230
	Ampicillin	AMP 10				
	Cefuroxime	CXM 30	✓	✓		
	Ceftazidime	CAZ 10				
	Trimethoprim-Sulfamethoxazole	SXT 25	✓	✓		
	Gentamicin	CN 10	✓	✓		
	Nitrofurantoin	F 100		✓		
	Norfloxacin	NOR 10	✓			
ENTEROBACTERIA 1	Amikacin	AK 30	✓	✓	100 Multodiscs	95220
	Aztreonam	ATM 30	✓	✓		
	Piperacillin-Tazobactam	TZP 36		✓		
	Chloramphenicol	C 30	✓	✓		
	Trimethoprim-Sulfamethoxazole	SXT 25	✓	✓		
	Netilmicin	NET 10		✓		
	Cefotaxime	CTX 5		✓		
	Nitrofurantoin	F 100		✓		
ENTEROBACTERIA 2	Levofloxacin	LEV 5	✓	✓	100 Multodiscs	95240
	Cefoxitin	FOX 30	✓	✓		
	Piperacillin-Tazobactam	TZP 36		✓		
	Cefuroxime	CXM 30	✓	✓		
	Streptomycin	S 300				
	Netilmicin	NET 10		✓		
	Ampicillin-Sulbactam	AMS 20	✓	✓		
	Tobramycin	TOB 10	✓	✓		



Multodisc with eight antimicrobial agents

*Customized
 antibiotics &
 concentrations
 welcome!*



page 2 of 2

Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.

Description	µg	CLSI	EUCAST	Packaging	Ref.	
PSEUDOMONAS	Amikacin Ticarcillin Piperacillin-Tazobactam Ceftazidime Gentamicin Netilmicin Levofloxacin Tobramycin	AK 30 TC 75 TZP 36 CAZ 10 CN 10 NET 10 LEV 5 TOB 10	✓ ✓ ✓ ✓ ✓ ✓ ✓ ✓	✓ ✓ ✓ ✓ ✓ ✓ ✓ ✓	100 Multodiscs	95250
STAPH	Cefoxitin Levofloxacin Netilmicin Tetracycline Erythromycin Gentamicin Rifampicin Linezolid	FOX 30 LEV 5 NET 10 TE 30 ERY 15 CN 10 RD 5 LNZ 10	✓ ✓ ✓ ✓ ✓ ✓ ✓ ✓	✓ ✓ ✓ ✓ ✓ ✓ ✓ ✓	100 Multodiscs	95260
STREPTO	Penicillin G Levofloxacin Teicoplanin Vancomycin Erythromycin Tetracycline Linezolid Rifampicin	P 10 IU LEV 5 TEC 30 VA 5 E 15 TE 30 LNZ 10 RD 5	✓ ✓ ✓ ✓ ✓ ✓ ✓ ✓	✓ ✓ ✓ ✓ ✓ ✓ ✓ ✓	100 Multodiscs	95290
ACINETOBACTER	Doripenem Meropenem Levofloxacin Gentamicin Trimethoprim-Sulfamethoxazole Tobramycin Ciprofloxacin Netilmicin	DOR 10 MRP 10 LEV 5 CN 10 SXT 25 TOB 10 CIP 5 NET 10	✓ ✓ ✓ ✓ ✓ ✓ ✓ ✓	✓ ✓ ✓ ✓ ✓ ✓ ✓ ✓	100 Multodiscs	95270
YEASTS	Caspofungin Fluconazole Posaconazole Voriconazole Amphotericin B Ketoconazole Flucytosine Nystatin	CAS 5 FLU 25 POS 5 VO 1 AMB 20 KCA 10 AFY 1 NY 100 IU	✓ ✓ ✓ ✓ ✓ ✓ ✓ ✓	✓ ✓ ✓ ✓ ✓ ✓ ✓ ✓	100 Multodiscs	95280

