

ANTIBIOTIC AGAR PENI CYLINDER DISPENSER



Stainless Steel Peni Cylinder Dispenser holds 300+ cylinders and can dispense 4 or 6 cylinders per plate. Cylinders are automatically and accurately dropped onto a Petri dish. Specify Petri plate size and type (glass or plastic) when ordering.

- *Dispensers manufactured from Stainless Steel and Aluminum*
- *Dispenser dimensions:
13 x 13 x 31 cm*
- *Cylinders manufactured from Stainless Steel*
- *Cylinder dimensions:
Inner diameter: 6 mm
Outer diameter: 8 mm
Length: 10 mm*
- *100 Stainless Steel Peni Cylinders per box.*



Giles Scientific Inc.
Santa Barbara, CA
Phone: 805-963-3876
Fax: 805-963-7768
E-mail: giles@biomic.com
www.biomic.com