

Biological Indicators

for monitoring and validating the sterilization process

Self-contained for steam sterilization

Product with minimum request



Description	Packaging	Ref.
STERILtest GST E6 self-contained Devices with nutrient medium, pH indicator and discs with spores of <i>G. stearothermophilus</i> (1-5x10 ⁶ CFU/disc).	100 units	91100
STERILtest GST E5 self-contained Devices with nutrient medium, pH indicator and discs with spores of <i>G. stearothermophilus</i> (1-5x10 ⁵ CFU/disc).	100 units	91101 #

Ampoules for steam sterilization



Description	Packaging	Ref.
STERIL CONTROL GST E6 ampoules Ampoules with nutrient medium, pH indicator and spores of <i>G. stearothermophilus</i> (1-5x10 ⁶ CFU/ampoule).	10 Ampoules 50 Ampoules	91010 91050
STERIL CONTROL GST E5 ampoules Ampoules with nutrient medium, pH indicator and spores of <i>G. stearothermophilus</i> (1-5x10 ⁵ CFU/ampoule).	10 Ampoules 50 Ampoules	91051 91052

Strips for steam sterilization

Product with minimum request



Description	Packaging	Ref.
STRIP CONTROL GST E6 Strips with spores of <i>G. stearothermophilus</i> (1-5x10 ⁶ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91055
STRIP CONTROL GST E5 Strips with spores of <i>G. stearothermophilus</i> (1-5x10 ⁵ CFU/strip) and tubes of nutrient broth with pH indicator.	100 strips 20 strips 20 tubes	91155 91056 #
STRIP CONTROL GST E4 Strips with spores of <i>G. stearothermophilus</i> (1-5x10 ⁴ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91057 #
STRIP CONTROL GST E3 Strips with spores of <i>G. stearothermophilus</i> (1-5x10 ³ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91058 #

Spores suspensions for steam sterilization

Description	Packaging	Ref.
SUSPENSION SPORE GST E8 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁸ CFU/mL).	10 mL vial	82002
SUSPENSION SPORE GST E7 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁷ CFU/mL).	10 mL vial	82003
SUSPENSION SPORE GST E6 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁶ CFU/mL).	10 mL vial	82004
SUSPENSION SPORE GST E5 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁵ CFU/mL).	10 mL vial	82005

Strips for dry heat or ethylene oxide sterilization

Product with minimum request

Description	Packaging	Ref.
STRIP CONTROL BAT E8 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁸ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91059 #
STRIP CONTROL BAT E7 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁷ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91062 #
STRIP CONTROL BAT E6 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁶ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91063
STRIP CONTROL BAT E5 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁵ CFU/strip) and tubes of nutrient broth with pH indicator.	100 strips 20 strips 20 tubes	91067 91064 #
STRIP CONTROL BAT E4 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ⁴ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91065 #
STRIP CONTROL BAT E3 Strips with spores of <i>B. atrophaeus</i> (1-5x10 ³ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91066 #

Spore suspension for dry heat or ethylene oxide sterilization



Description	Packaging	Ref.
SUSPENSION SPORE BAT E10 Sterile suspension of <i>B. atrophaeus</i> spores in distilled water (1-5x10 ¹⁰ CFU/mL).	10 mL vial	82006
SUSPENSION SPORE BAT E9 Sterile suspension of <i>B. atrophaeus</i> spores in distilled water (1-5x10 ⁹ CFU/mL).	10 mL vial	82007
SUSPENSION SPORE BAT E8 Sterile suspension of <i>B. atrophaeus</i> spores in distilled water (1-5x10 ⁸ CFU/mL).	10 mL vial	82008
SUSPENSION SPORE BAT E7 Sterile suspension of <i>B. atrophaeus</i> spores in distilled water (1-5x10 ⁷ CFU/mL).	10 mL vial	82009
SUSPENSION SPORE BAT E6 Sterile suspension of <i>B. atrophaeus</i> spores in distilled water (1-5x10 ⁶ CFU/mL).	10 mL vial	82010
SUSPENSION SPORE BAT E5 Sterile suspension of <i>B. atrophaeus</i> spores in distilled water (1-5x10 ⁵ CFU/mL).	10 mL vial	82011

Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.

Self-contained for ethylene oxide sterilization

Product with minimum request



Description	Packaging	Ref.
STERILtest BAT E6 self-contained Devices with nutrient medium, pH indicator and discs with spores of <i>B. atrophaeus</i> (1-5x10 ⁶ CFU/disc).	100 units	91150
STERILtest BAT E5 self-contained Devices with nutrient medium, pH indicator and discs with spores of <i>B. atrophaeus</i> (1-5x10 ⁵ CFU/disc).	100 units	91151 #

Crusher for self-contained

Description	Packaging	Ref.
CRUSHER Steel pincer for the easy break of self-contained vials.	1 pincer	96790

Strips for irradiation sterilization

Product with minimum request

Description	Packaging	Ref.
STRIP CONTROL PUM E7 Strips with spores of <i>B. pumilus</i> (1-5x10 ⁷ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91071 #
STRIP CONTROL PUM E6 Strips with spores of <i>B. pumilus</i> (1-5x10 ⁶ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91072
STRIP CONTROL PUM E5 Strips with spores of <i>B. pumilus</i> (1-5x10 ⁵ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91073 #
STRIP CONTROL PUM E4 Strips with spores of <i>B. pumilus</i> (1-5x10 ⁴ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91074 #
STRIP CONTROL PUM E3 Strips with spores of <i>B. pumilus</i> (1-5x10 ³ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91075 #

Steel coupons for hydrogen peroxide sterilization

Product with minimum request



Description	Packaging	Ref.
OXI CONTROL E6 SS Steel coupons with spores of <i>G. stearothermophilus</i> (1-5x10 ⁶ CFU/coupon) and tubes of nutrient broth with pH indicator.	20 coupons 20 tubes	91082
OXI CONTROL E5 SS Steel coupons with spores of <i>G. stearothermophilus</i> (1-5x10 ⁵ CFU/coupon) and tubes of nutrient broth with pH indicator.	20 coupons 20 tubes	91083 #

Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.

Poliflex coupons for hydrogen peroxide sterilization

Product with minimum request



Description	Packaging	Ref.
OXI CONTROL E6 PLX Poliflex coupons with spores of <i>G.stearothermophilus</i> (1-5x10 ⁶ CFU/coupon) and tubes of nutrient broth with pH indicator.	20 coupons 20 tubes	91085
OXI CONTROL E5 PLX Poliflex coupons with spores of <i>G.stearothermophilus</i> (1-5x10 ⁵ CFU/coupon) and tubes of nutrient broth with pH indicator.	20 coupons 20 tubes	91086 #

Strips for hydrogen peroxide sterilization

Product with minimum request

Description	Packaging	Ref.
OXI CONTROL E6 STRIP Strips with spores of <i>G.stearothermophilus</i> (1-5x10 ⁶ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91088
OXI CONTROL E5 STRIP Strips with spores of <i>G.stearothermophilus</i> (1-5x10 ⁵ CFU/strip) and tubes of nutrient broth with pH indicator.	20 strips 20 tubes	91089 #

Spore suspensions for hydrogen peroxide sterilization

Description	Packaging	Ref.
SUSPENSION SPORE GST VHP E8 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁸ CFU/mL).	10 mL vial	82030
SUSPENSION SPORE GST VHP E7 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁷ CFU/mL).	10 mL vial	82031
SUSPENSION SPORE GST VHP E6 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁶ CFU/mL).	10 mL vial	82032
SUSPENSION SPORE GST VHP E5 Sterile suspension of <i>G. stearothermophilus</i> spores in distilled water (1-5x10 ⁵ CFU/mL).	10 mL vial	82033

PVC coupon for hydrogen peroxide sterilization

Description	Packaging	Ref.
OXI CONTROL E6 PVC PVC coupons with spores of <i>G.stearothermophilus</i> (1-5x10 ⁶ CFU/coupon) and tubes of nutrient broth with pH indicator.	20 coupons 20 tubes	91092

Aluminum coupons for hydrogen peroxide sterilization

Description	Packaging	Ref.
OXI CONTROL E6 ALUMINUM Aluminum coupons with spores of <i>G.stearothermophilus</i> (1-5x10 ⁶ CFU/coupon) and tubes of nutrient broth with pH indicator.	20 coupons 20 tubes	91098

Glass coupons for hydrogen peroxide sterilization

Description	Packaging	Ref.
OXI CONTROL E6 GLASS Glass coupons with spores of <i>G.stearothermophilus</i> (1-5x10 ⁶ CFU/coupon) and tubes of nutrient broth with pH indicator.	20 coupons 20 tubes	91096

Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.



Via Scozia zona ind.le, 64026 Roseto degli Abruzzi (Te) Italy

Tel. +39 0858930745 Fax +39 0858930330 www.liofilchem.net

liofilchem@liofilchem.net

ref. 6551066
Rev. 3 - 04.03.2013
© Liofilchem 2013
Printed in Italy